

**TECHNICAL BULLETIN: TESTING COUNT AND POST-VAS PROFICIENCY SAMPLES ON SQA SYSTEMS**  
**SQA-VISION / SQA-V ALL VERSIONS | July 22, 2020**

**GENERAL OVERVIEW:**

1. Allow specimens, and any reagents used, to warm to room temperature (at least 30 minutes) before testing. Assume this is a liquefied specimen.
2. To dislodge the pellet, **do not vortex**. Aspirate the sample up and down 5 - 6 times using a transfer pipette (taking care to keep the pipette below the surface of the sample).
3. Run the analyzer to perform a count per the instructions below (based on instrument type). If no sperm are present, report 0.0 M/mL.
4. For **Fertility** specimens, indicate whether the obtained sperm count would be considered normal or abnormal based on your laboratory's reference range.

**TESTING INSTRUCTIONS:**

**SQA-V Users:**

1. **Post Vasectomy** samples should be run on the standard Post Vas mode of the SQA-V system per user guide instructions. The optional "Manual Capillary Scan" will need to be performed per the SQA-V's instructions to ensure an accurate proficiency test.
2. **Fertility** samples should be run on QC (Stabilized Sperm) mode:
  - a. Set up the SQA-V Controls default to Stabilized Sperm and apply in V-Sperm: **SET-UP > SQA-V > SQA-V DEFAULTS > Select "Stabilized Sperm" > Enter "25" as the Target > Enter "25" As the Range > Apply.**
  - b. Load the SQA-V Testing capillary and run the specimen as a control. On the SQA-V/SPERMALITE, select: **RUN CONTROLS > LEVEL 1.**
  - c. *Note: Do not run specimen as a patient or latex bead QC sample.*

If the stabilized sperm sample run on the SQA-V GOLD in the stabilized sperm Quality Control mode results in ZERO (below the 2 M/mL dynamic range), re-run the sample in the Fresh mode as outlined below:

- a. Turn on the SQA-V GOLD system and wait until auto-calibration/self-testing is completed.
- b. Go to: MAIN MENU>TEST NEW PATIENT and enter:  
PATIENT ID: **Sample #**  
BIRTH DATE: **Test date**
- c. From the next screen select:  
SAMPLE TYPE: **FRESH**  
WBC CONC: **< 1M/mL**



- d. Select “Yes” when asked “IS SAMPLE VOLUME SUFFICIENT FOR COMPLETE RETESTING  $\geq$  0.5 mL?”
- e. Repeat testing.

### SQA-Vision Users:

**Post Vasectomy** samples should be run on the standard Post Vas mode of the SQA-Vision system per user guide instructions. The “Manual Slide Scan” PV test type should be performed to ensure accurate results. Do not report automated results only.

**Fertility** samples should be run on Proficiency Mode:

#### Setting up the defaults:

- a. From the MAIN MENU of the VISION PC, select SETTINGS > PROFICIENCY.
- b. Enter a SAMPLE ID, DATE, and a NOTE (if necessary).
- c. Click SAVE.

#### Testing proficiency samples on the SQA-Vision:

- a. Select QC/PROFICIENCY from the main MENU and touch TEST NOW on the desired option.
- b. Mix the sample well by aspirating it in and out 10 times with a transfer pipette. **Avoid bubbles.**
- d. Load the SQA-V testing capillary and check **closely** for bubbles then wipe the capillary tip clean.
- c. Insert capillary into the testing chamber when prompted to do so.
- e. Report the results received.
- d. Repeat the 5 steps above for all samples.

CONTACT MES, LLC FOR ADDITIONAL INFORMATION AND QUESTIONS: [service@mes-llc.com](mailto:service@mes-llc.com).

Distribution: All SQA-V and Vision Users

Authority: Eric Carver, General Manager

*“Remember, it ALL Started with a Sperm!”*



[www.mes-global.com](http://www.mes-global.com)