

Technical Release Bulletin

SQA-V Capillary Blue Valves (pistons)

Issue date: Monday, August 13th, 2007
Bulletin number: 13_AUG_2007_SQAV Capillary_Piston
Subject: Defective Blue Valves in the SQA-V Testing Capillary

Status:

There is reported that the Blue Valve (piston) of the SQA-V capillaries are difficult to close.

Troubleshooting:

- Advise the laboratory to push VERY hard on the blue valve prior to running a test in order to completely and effectively close the valve.



Explanation/Possible Causes:

- A new vendor supplied the dark blue valves (pistons).
- The rubber used in the dark blue pistons is softer and more porous than the previous rubber (light blue).
- The more porous and softer rubber absorbs the silicone sprayed on the pistons more completely causing more difficulty and more pressure required to close the piston.

Recommendations:

- Replace the dark blue softer rubber pistons with new pistons supplied by the manufacturer.

Instructions for Replacement of the dark blue valves (pistons):

<p>Remove all defective dark blue valves (pistons) from the existing lots of finished capillaries as advised by MES.</p> <p>Note the structure of the blue valve (piston)</p> <ul style="list-style-type: none"> • Rib _____ • Hole _____ 	
<ol style="list-style-type: none"> 1. Spray special silicon (for plastic and rubber)- on the blue valves (pistons) This has been done by MES and we have provided EXTRA silicon spray in case you notice the pistons are not going into the hole easily. 2. Inset the blue valve (piston) into open hole on the top side of the Capillary Runner as follows: <ul style="list-style-type: none"> • Align the rib of the blue valve (piston) with the slot in the section of the Capillary Runner where the valve is inserted. • Line-up the holes of the blue valve (piston) with the holes in the Capillary Runner. • Push the blue valve (piston) in all the way. 	
<ol style="list-style-type: none"> 3. Using the jig, push the blue valve (piston) to its proper position in the capillary: <ul style="list-style-type: none"> • Place the capillary on the top of the jig with the piston hole on the top of the point of the jig • Press the point of the jig into the capillary hole and move the blue valve (piston) into the correct position. 4. The Capillary Separating Valve will be about 1 mm above the top of the Capillary Runner, with the holes in alignment with the Capillary Runner when finished. 	