



## Hardware and Software Requirements for SQA-V (all systems) and X-Sperm

The following tables define the requirements and settings for optimal SQA-V and X-Sperm performance:

Computer	Desktop and Laptop computers have the same requirements
CPU	Intel Core I5 CPU 2.4GHz Or better
RAM	4GB
Video Card	3D Video Cards to support high resolution (16:10 – 1440 x 900)
Screen resolution	Wide screen resolution. The screen should support resolution of 16:10 or 16:9 (i.e. 1440 x 900)
USB Ports	<p><b>Four free</b> available USB ports</p> <p>First USB port is for USB to RS232 converter to communication between the PC and the SQA</p> <p>Second USB port is for USB to RS232 converter to sending data over the LIS system</p> <p>Third USB port is for the security key (X-Sperm)</p> <p>Fourth port is for connecting the video capture device</p>
CD-ROM	Standard CD-ROM compatible drive
Hard Disk	At least 300 GB free hard drive for capturing videos and images
Keypad	Standard US 101 Keypad
Mouse	Standard mouse
Components	
Video Capture device + Driver	Supplied by the manufacturer with each X-Sperm software package
RCA video cable	Supplied by the manufacturer with each X-Sperm software package
Security key	Supplied by the manufacturer with each X-Sperm software package
TWO RS232 cable	Two Standard RS232 cables: One for connecting the SQA-V and one to connect to the LIS host



Software	
X-Sperm	Installation disk Supplied by the manufacturer
Office 2003 or higher	Some software controls/functions require Office (Windows 7) – WORD and EXCEL required for downloading test results.
Operating system	Windows 7, Windows XP
Video Capture Device Driver	Supplied by manufacturer

Settings	
The PC printer must be defined	In order to enable a report in the X-Sperm software, a default printer must be defined in the PC
Windows Classic style	The windows style must be set to "Classic"
User account control setting (UAC)	User account control setting must be set to lower security setting (Windows 7 32 bit only)
Display resolution (wide)	Set to 1440 x 900 (wide)
Visualization default settings - 1	The video settings should be adjust to Analog Video Standard - PAL_B Video subtype – default Video size – 640 x 480 Video input – composite
Visualization default settings - 2	Update the visualization control settings in order to get a good image
Video compression tool	Set to Div-X MPEG 4

	Name	Title	Signature	Date
Written	Beni Cohen	Director of R&D		4-MAY-14
Approved	Marcia Deutsch	CEO		4-MAY-14
Approved	Dr. Lev Rabinovitch	CTO		4-MAY-14