BSE B-Sperm™ & SQA-Vb

User Guide

Catalog # 7410 Version 2.0 May 2015



Section	TABLE OF CONTENTS	Page
1	Overview	 3
2	SQA-Vb Overview	 3
3	Start-up: SQA-Vb	 5
4	Start-up: B-Sperm	 5
5	B-Sperm Navigation Overview	 5
6	Test Bull	 6
7	Test Data	 10
	BSE Results	 12
	QwikClick [™] Morphololgy Counter	 12
	Test Report/Morphology Report	 15
	Data Analysis Features	 16
8	Video	 16
9	Controls	 17
10	Export	 19
11	Set-up	 19
	BSE Report	 19
	QwikClick™ Morphology Counter Set-up	 20
	Defaults: Cut-offs/Normal ranges/Failure codes	 20
	Language	 21
	Password	 22
	Port	 22
	SQA-Vb	 22
12	EXIT	 23
Appendix	1: Semen Sample Preparation	 24
	2: Capillary Filling Instructions	 25
	3: SQA-Vb Cleaning Instructions	 26
	4: Capillary Washing/Drying Instructions	 27
	5: The Visualization System	 28
	6: Sample Heating Options	 29
	7: Glossary of Terms	 30
	8: SQA-Vb System Specifications	 31
	9: Product Performance Data	 34
	10: SQA-Vb Service Menu	 38
	Add I-button tests to counter	 38
	Date Setting	 38
	11: Troubleshooting: SQA-Vb Warning Screens	 39

Section 1: Overview

The SQA-Vb and B-Sperm[™] bull BSE management software work together to provide an integrated system for bull breeding soundness evaluation (BSE).

Both the **physical examination and semen testing** elements of the BSE are included in the software and integrated into a final BSE report.

FRESH bull semen samples are prepared and aspirated into a multi-use, washable testing capillary. The capillary is inserted into the measurement compartment of the SQA-Vb where state of-the-art technology in electro-optics, computer algorithms and video microscopy produce rapid, accurate and precise analysis of the sample. Samples can be viewed using the on-board video visualization system with a magnification range of X300 through X500.

All bull/sample data entry, settings and reports are handled in the PC-based BSE management software. The following functions and features of the system are listed below:

- Rapid data entry that is updated during the testing process or from pre-set tables.
- Fast test results Less than one minute is required to report all semen parameters. Perfect for high throughput.
- Easily assess morphology defects using the QwikClick[™] feature.
- Herd data analysis (histograms).
- Numeric and graphical Bull information/reports.
- Bull semen samples can be viewed on a large PC screen (video or image) and images can be captured and stored with individual bull records.
- Report options: BSE, Test Report, Morphology (differential).

Section 2: SQA-Vb Overview







Sect	Before tes settings in Turn on b	sting for the find the second strength for the find the second se	irst time or i lease refer to Vb and B-Sp	n the event o the Set-up erm to begin	of a change, section). BSE testing	update the de
Click t	he Test B	ull button and	I the screen	below will be	displayed:	
	B-Sperm Manageme	ent System				
	Advanted April Provide the Provide State	B-Sp	erm Version 182.0	11		
		Test Bull				
	Test Bull	Owner Name	Address	Telephone	Fax	E - mail
	Test Data Video	MES .	20 Alon Hatavor	04-6373981	04-6373984	mes@mes-ltd.con
	Controls	Bull ID v	Bull Name	Birth Date	Age(Months)	Breed
	Export Set-Up	*				
	Exit	123	Star	01/04/2012	13 Months	Holstein
	l					
	l	Semen Examination	Physical Examination	Delete	Print]
EN	<u>ter</u> :	Semen Examination	Physical Examination	Delete	Print]
<u>EN</u> •	<u>TER</u> : OWNER	Semen Examination	Physical Examination	Delete	Print D a given ow	ner (first row)
<u>EN</u> • •	TER: OWNER BULL ID,	Semen Examination NAME and oth , NAME (if rec	Physical Examination her informat quired), BIR	Delete ion related to TH DATE (A	Print o a given ow GE will be ca	ner (first row) alculated
EN • • The	TER: OWNER BULL ID, automatic	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be	Physical Examination her informat quired), BIR EED.	Delete ion related to TH DATE (A each bull is t	Print o a given ow GE will be constructed OR in	ner (first row) alculated advance by e
EN • The all 1	TER: OWNER BULL ID, automatio e bull infor the bulls ir	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table.	Physical Examination her informat quired), BIR EED. e entered as	Delete ion related to TH DATE (A each bull is t	Print o a given ow GE will be content cested OR in	ner (first row) alculated advance by e
EN • The all t	TER: OWNER BULL ID, automatio e bull infor the bulls in LETE:	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table.	Physical Examination her informat quired), BIR E ED . e entered as	Delete ion related to TH DATE (A each bull is t	Print o a given ow GE will be co cested OR in	ner (first row) alculated advance by e
EN • The all f DE	TER: OWNER BULL ID, automatic bull infor the bulls in LETE: Bull record	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table. ds by selecting putton	Physical Examination her informat quired), BIR EED. e entered as g the record	Delete ion related to TH DATE (A each bull is t in the far lef	Print o a given ow GE will be co cested OR in t hand box a	ner (first row) alculated advance by e and pressing t
EN • • The all 1 • •	TER: OWNER BULL ID, automatic bull infor the bulls in LETE: Bull record DELETE b Owner by	Semen Examination NAME and oth , NAME (if red cally) and BRI mation can be n this table. ds by selecting button.	Physical Examination her informat quired), BIR EED. e entered as g the record small square	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ov	Print o a given ow GE will be co cested OR in it hand box a wner Name a	ner (first row) alculated advance by e and pressing t
EN • The all t DE •	TER: OWNER BULL ID, automatic e bull infor the bulls in LETE: Bull recor DELETE to Owner by DELETE to	Semen Examination NAME and oth , NAME (if rec cally) and BRI rmation can be n this table. rds by selecting outton. r selecting the poutton. All the	Physical Examination her informat quired), BIR ED. e entered as g the record small square e bull recor	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ov ds related t	Print D a given ow GE will be co cested OR in t hand box a wner Name a o this owne	ner (first row) alculated advance by e and pressing t and pressing t er will be de
EN • • • • • Pre	TER: OWNER BULL ID, automatic bull infor the bulls in LETE: Bull recor DELETE to Owner by DELETE to ss the PR:	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be mation can be mation can be n this table. ds by selecting button. selecting the button. All the INT button to	Physical Examination her informat quired), BIR E D . e entered as g the record small square bull recor print the bu	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ou ds related t	Print D a given ow GE will be co cested OR in it hand box a wner Name a o this own le for a give	ner (first row) alculated advance by e and pressing t and pressing t er will be del n Owner.
EN • • • • • Pre If th of t	TER: OWNER BULL ID, automatic e bull infor the bulls in LETE: Bull record DELETE to Owner by DELETE to ss the PR: he OWNE the three of	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table. This table.	Physical Examination her informat quired), BIR ED . e entered as g the record small square bull recor print the bun there twice lained in the	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ou ds related t ill setting tab e, the message message.	Print o a given ow GE will be co sested OR in it hand box a wher Name a o this owno le for a give ge below is c	ner (first row) alculated advance by e and pressing t and pressing t er will be del n Owner. displayed. Sel
EN • • • • Pre If the of the •	TER: OWNER BULL ID, automatic bull infor the bulls in LETE: Bull record DELETE to Somer by DELETE to so the PR: the OWNE the three co	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table. This table. T	Physical Examination her informat quired), BIR ED. e entered as g the record small square bull recor print the bun there twice lained in the	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ou ds related t ill setting tab e, the message	Print o a given ow GE will be co sested OR in it hand box a wher Name a o this owne le for a give ge below is c	ner (first row) alculated advance by e and pressing t and pressing t er will be del n Owner. displayed. Sel
EN • • • • Pre If tl of t	TER: OWNER BULL ID, automatic e bull infor the bulls in LETE: Bull record DELETE to Owner by DELETE to ss the PR: he OWNE the three of	Semen Examination NAME and oth , NAME (if rec cally) and BRI mation can be n this table. This table. This by selecting outton. Selecting the button. All the INT button to R (Name) is in options as exp This OWNER In order to p o VES - b o NO-to o NO-to	Physical Examination her informat quired), BIR ED. e entered as g the record small square print the built n there twice lained in the already exists! Do you way occeed, click: preade a new owner with update the existing own- - the optimum authout for a	Delete ion related to TH DATE (A each bull is t in the far lef e near the Ou ds related t ill setting tab e, the message message.	Print D a given ow GE will be ca sested OR in it hand box a wher Name a o this owne le for a give ge below is c with the same name?	ner (first row alculated advance by e and pressing t and pressing t er will be de n Owner. displayed. Se

START TESTING:

Select the desired **BULL** and click on **Semen Examination** or **Physical Examination** to start the testing process. •

SEMEN EXAMINATION: Semen Click **SEMEN EXAMINATION** to display the screen below: Examination B-Sperm Manage _ 8 × ATELIS **B-Sperm** rsion 182.011 Test Bull > Semen Examination Test Bull Examination Date Test Data 01/05/2013 Video Controls Owner Bull ID Bull Name Breed Birth Date Age (Months) ort MES Sta Holstein 01/04/2012 13 Months Semen Appearance Normal • Semen Volume (ml) 1 Sample # Physical Cancel Continue >> Examination Select/Enter: SEMEN APPEARANCE: Normal or Semi-Transparent (Thin). SEMEN VOLUME (ml) SAMPLE # Click **CONTINUE>>** to begin testing or **CANCEL** to cancel the test. Click **PHYSICAL EXAMINATION** to open the input screen for the BSE report. FRESH SAMPLE PREPARATION AND TESTING: NOTE: When If **SEMEN APPEARANCE: Normal** is selected, after clicking **Continue>>**, the preparing following sample preparation screens are displayed: FRESH samples with Prepare a Fresh sample for testing: Prepare a Fresh sample for testing: Normal Appearance 1. Diluent: 2 ml (pre-heated to 37⁰C) 1. Diluent: 2 ml (pre-heated to 37⁰C) for testing the 2. Semen: 200 ul 2. Semen: 200 ul diluent volume 3. Mix sample thoroughly and filter impurities 3. Mix sample thoroughly and filter impurities is alwavs 2 4. Fill and wipe a pre-heated testing capillary 4. Fill and wipe a pre-heated testing capillary milliliters. Do not touch the SQA-Vb at this time **Ready for testing!** Insert the testing capillary into the SQA-Vb Autocalibration is in progress NOTE: DO NOT Prepare a FRESH sample based on the instructions. pre-heat the Prepare 2ml aliquots of diluent pre-heated to 37^oC (98.6^oF). FRESH 0 sample! Only Pre-heat the testing capillaries to $37^{\circ}C$ (98.6°F). 0 pre-heat the diluent and DO NOT HEAT THE SAMPLE! ONLY THE DILUENT AND CAPILLARY or 0 the testing motility will be impacted! capillaries or motility will be If SEMEN APPEARANCE: Semi-Transparent (Thin) is selected and the SEMEN impacted. **VOLUME => 1 ml**, after clicking **Continue>>** the sample preparation screens below are displayed: NOTE: Pre-heat the Prepare a Fresh sample for testing: Prepare a Fresh sample for testing: FRESH 1. Filter impurities 1. Filter impurities sample per 2. Heat sample for 1 minute 2. Heat sample for 1 minute onscreen 3. Fill and wipe a pre-heated testing capillary 3. Fill and wipe a pre-heated testing capillary instructions when run without Do not touch the SQA-Vb at this time Ready for testing!

Insert the testing capillary into the SQA-Vb

Autocalibration is in progress

dilution.



- The results will be saved by the system.
- The QwikClick Morphology Counter can automatically be displayed (if this is preset). Alternatively it can be opened offline for each bull record when needed.



Click the box next to **View entire history** to view past BSE records for the bull.

- The History section data shows if the bull was classified as **Satisfactory or Unsatisfactory Potential Breeder** or classification was **Deferred**.
- Click: **SAVE** to store the results in the BSE archive.

Section 7: Test Data

Test Data

Click **Test Data** option from the main menu and the following grid will be displayed:

ATE		stem	в	-Spern		ersion 182.0	.1						
	Test Dat BSE	ta Results	Qv	ikclick	Test Report	t Mor	ph. Report	Dat	ta Analysis	← BACK			
Test Bull	From D	Date 2 /1	9/2013	_ To Date	2 /19/2013		Colorat	-					
Test Data	Numbe	r of Reco	ords 3		10/2010	<u> </u>	Select				Sort	Hide	View All
Video Controls		Test	_	Bull	Bull		_	Case	Sample	Semen	Sperm	Motility	Prog.
Video Controls Export		Test Date [▽]	Time ⊽	Bull Name	Bull ID	Breed	Owner	Case #	Sample #	Semen Volume [ml]	Sperm Conc. [M/ml]	Motility [%]	Prog. Motility [%]
Video Controls Export Set-Up	▶ 2/	Test Date [▽] 19/2013	Time ⊽ 10:50	Bull Name Ace	Bull ID 125	Breed Holstein	Owner John	Case #	Sample #	Semen Volume [ml] 3.500	Sperm Conc. [M/ml] 338.2	Motility [%] 79.0	Prog. Motility [%] 77.0
Video Controls Export Set-Up	▶ 2/ 2/	Test Date ∇ 19/2013 19/2013	Time ⊽ 10:50 10:43	Bull Name Ace King	Bull ID 125 124	Breed Holstein Holstein	Owner John John	Case #	Sample # 1	Semen Volume [ml] 3.500 4.500	Sperm Conc. [M/ml] 338.2 350.0	Motility [%] 79.0 82.3	Prog. Motility [%] 77.0 80.5

Click in the grey box of the first column to select one record or drag the arrow \blacktriangleright to select multiple records. The following information is displayed in the **Test Data** table:

Sample/Bull Information:

- **Date** The date the sample was tested.
- **Time** The time the sample was tested.
- Bull Name Entered through the Test Bull data entry screen (not required).
- **Bull ID** A unique bull identification # designated by the owner (required).
- **Breed** The bull breed entered in the data entry screen.
- **Owner** The name of the bull owner.
- **Case #** A number designated by the user.
- Sample # The sample number designated by the user.
- Semen Volume The ejaculate volume (ml).

<u>Test Results</u>:

- **Semen parameters**: Sperm Concentration, Motility, Progressive Motility, % Normal Cells, % Primary Abnormalities and % Secondary Abnormalities.
- Sample Information: Sample Pass/Fail. If all parameters pass the pre-set cutoff, Sample Pass is indicated. The Sample Fail will be displayed if at least one parameter is below the cut-off. The Sample Fail can be changed to Sample Pass manually if desired.

Additional Information:

- A variety of the following icons indicate different features:
- **Graphs**: This icon will appear when there are enough test results to graph. Click on the icon and select the parameter desired to produce a graph.
- <> Test results out of normal range: This icon alerts the user that at least one test result is **out of the normal range.** Refer to the Set-up>Normal Ranges section for details).
- Picture indicator: This indicates that images are attached to a record. Click to view images captured from stained morphology slides (can be attached to each test).

Video Indicator: This Click to view. Only on	icon indicates a video clip has been attached to the record. e video clip can be attached to each test record.
Three buttons appear at the	e bottom of the Test Data screen:
Clea	ar All Select All Delete
• Clear All: Click to des	select all selected records.
Select All: Click to se	ect all records.
Delete: Click to delet	e selected record.
Three buttons are located in	n the upper right hand corner of the Test Data screen.
	Sort Hide View All
• Click: SORT and then sort (Date, Bull ID, Sp	click on the column header that contains the information to perm Conc., etc.).
• Click: HIDE and then	select the columns to hide. The columns will disappear.
Click: VIEW ALL to re	e-activate ALL of the hidden columns.
Five buttons and a date ran screen: Test Data BSE Results Qwike	Ige setting are displayed at the top of the Test Data
From Date 2 /18/2013	▼ To Date 2 /18/2013 ▼ Select
Set the date range to display however this can be re-set by following message is shown:	¹ data within that timeframe. The default is one month, y the user. If the selected date range is too broad, the
P-S	Sperm
Т	he date range selected will take several minutes to retrieve. You'd you like to releat a different date range?
v	Yes No

Te	st Data > BSE Results			Semen Exa	mination
	BSE Results 🗲 BACK			Collection Method	EE 💌
	BULL BREEDING SOU	NDNESS EVALUATION]	Response	Ejaculation 💌
0	wner	John		Semen Char	acteristics
A	ddress	20 Goldberg st, LA		Concentration (M/ml)	986.4
Т	elephone	98735612		Motility (%)	43.6
F	ах	98735613		Progressive Motility (%)	42.5
E	-Mail	john@gmail.com		% Normal Cells	87.4
C	ase #	1		% Primary Abnormalities	4.7
E	xamination Date	2/17/2013		WBC	4.7
B	ull Name	Star		BBC	
В	ull ID	123	-	Other	
	reed	Holstein		Classifi	cation
B	Inth Date	1/1/2012 12 Maatha	-	Interpretation of data resulting from	this examination would indicate
É É	ge (monuns) History : Pr	TO MONUNS		that <u>on this date</u> , this bull is a:	this examination would indicate
	ate	evious Dat			unden wie of w
	aic ase #			Satisfactory potential bre Re-examination recommended on	eder <u>v</u> v <u>v</u>
	lassification		-		
F	Physical E	kamination		I his buil has been examined to quality of semen only. Unless o	or pnysical soundness and therwise noted no
в	ody Condition Score	6 🔻		diagnostic tests were undertake	en for libido, mating ability or
в	ody Condition	Good T		infectious desease status of this	s bull.
P	elvic Height	12.8		Remarks and Interpretation (diagnos	sis, prognosis, recommendations
P	elvic Width	11.9			
P	elvic Area	152			
F	eet/Legs	 ✓ 			
E	yes	✓		Signed:	Uri
V	esicular Glands	 ✓ 	il l		MEMBER - SOCIETY FOR
A	mpullae/Prostate	 ✓ 			THERIOGENOLOGY
In	guinal Rings	 ✓ 		Clinic:	Medisoos
P	enis/Prepuce	 ✓ 			
Т	estes/Spermatic Cord	✓ •		Save	Report
E	pididymides	✓ •			
S	crotum (Shape)	 ✓ 			
0	ther				
5	crotal Circumference	30	_		
	Save	Beport			
		Hoport			
	Move the scr	lling button on th	o right	side of the arid to vi	iow all of the
•	information	ning bacton on a	ic right		
	iniornation.				
•	Owner and bu	Ill data, semen ar	nd physi	ical examination resu	ults are displayed
	BSE arid. Add	litional data such	as the h	oull classification, con	mments, etc. car
	manually ent	ared in the arid			
		ereu in the griu.			
•	Click: SAVE t	o retain the BSE	results.		
•	Click: Report	for the printable	BSE Re	eport.	
ow	KCLICK:				
			- -		
ihis	section describe	s now to use the	B-Sperr	m QwikClick' [™] Morph	ology Counter.
	INITIONS. The	owikClick™ foot	ure hac	the following option	s for accoccing h
			ure lids		s for assessing D
norr	phology if the fu	ily automated opt	ion is n	ot preferable:	
101	OwikClick™	Morpholoav Cou	inter (a	on-line right after the	automated test
101	z				
101 F			- /	tt line accoccy Tect	DataSOwikCliv
((• QwikClick™	Morphology Col	inter (o	on-line, access. Test	



The sample will **PASS** if all of the semen parameters have passed the cut-offs. The sample will **FAIL** if at least one semen parameter fails. The user can manually change the automatic PASS/FAIL indication.

0.0	Pass
30.0	Pass
25.0	Pass
70.0	Pass
	30.0 25.0 70.0

Automated Test Results				
Parameter	Results	Cut-off	Pass/Fail	
Sperm Conc. [M/ml]	1015.7	0.0	Pass	
Motility [%]	29.7	30.0	Fail	
Prog. Motility [%]	28.8	25.0	Pass	
Normal Cells [%]	61.9	70.0	Fail	
Sample Pass Reason for Failure	Fail		• ?	
Semen result	s) do not	pass cu	t-off!	

wikClick	k Morphology Counter			To use t
Keys	Classification	# Cells	7.	icor
F1	Normal Cells	0	0	Note
	Primary Abnormalities			Note.
F2	Underdeveloped	0	0	The fun
F3	Double Forms	0	0	defects
F4	Heads	0	0	set by t
F5	Midpiece	0	0	Proced
F6	Proximal Droplets	0	0	
F7	Tails	0	0	• Prep
	Secondary Abnormaliti	es		place
F8	Heads	0	0	com
F9	Abaxial Implantation	0	0	• Use
F10	Distal Droplets	0	0	cour
F11	Tails	0	0	abno
F12	Other	0	0	sper
21	Total	0	0.0	Sper
Sta	t Final Calculation	Field N	ew # 0	Click view SQA
				Click

e QwikClick Morphology Counter click

to display instructions.

tion keys that correspond to morphological re set by the manufacturer, but can be ree user in the B-Sperm set-up.

re:

- re a slide of non-motile sperm cells and it into the SQA-Vb visualization artment.
- ne PC function keys (F1, F2, etc.) for ing normal cells and morphological malities (one abnormality per natozoa).
- Next Field and move to another field of by slightly turning the silver knob of the Vb slide holder.
- FINAL CALCULATION to save the results.
- A differential Morphology assessment can be performed using stained smears or captured images.
- If the user exits this screen without counting, the **% Normal Cells** will be calculated based on automated results.
- If a differential Morphology is performed and the **Final Calculation** button is pressed, the % Normal Cells will be calculated based on the manual assessment.
- Differential Morphology results will be seen in the BSE and Morphology Report.

Capture Picture or Video. To save a picture or video clip and attach it to a record:



- Insert a slide with a semen sample into the visualization system of the SQA-Vb.
- Highlight a record in the Test Data grid to attach the video/picture (not necessary if the QwikClick[™] is opened on-line).



- Click: Stop Capturing to end the video capture process.
- Several pictures and one video can be saved and attached to the bull record.
- To **delete** a picture or video, select the desired image at the bottom table of the QwikClick screen (shown above) and click the **Delete** button.
- A picture or video can be deleted in the **Test Data** screen by clicking the picture or video icon in the **Test Data** grid and then clicking the Delete button (located at the bottom of the grid or video clip).



TEST REPORT:



Click on this icon to print a copy of the report.



Click on this icon to export the semen analysis report to Excel file. Click **Test Report** to view and print the following report displaying the bull data and semen analysis results for the selected period of time:



MORPHOLOGY REPORT:



Click **Morphology Report** to view and print the Morphology report, which displays the Sample Data and QwikClick differential morphology results of the selected records:

Mr	orpholog	y Report																			
	e	Zoor	n 100%	•																	
										SQA-\	/b										
							Ν	lor	oho	loa	v Repo	ort									
							R	eport	Date:	2/18/2	013 2:59:11										
Γ				Sar	mple Data					Normal Cells		Primo	ry Abno	rmalitics				Secondary	/ Abnorma	lities	
	Date	Time	Bull Name	Bull ID	Breed	Owner	Case #	Sampie #	Semen Volum c [m]	[%]	Underdeveloped [%]	Double Forms [%]	Heads [%]	Midpiece [%]	Proximal Droplets [%]	Tails [%]	Heads [%]	Abaxial Implantation [%]	Distal Droplets [%]	Tails [%]	Other [%]
$\left \right $	2/17/2013	16:50	Star	123	Holstein	.lohr	1	1	5.000	87 4	21	0.8	16	08	0.8	16	16	0.8	0.8	0.8	68





Controls Testing	Controls > Set-up
QwikCheck™	Level 1 Level 2 Negative Control
beads are	Lot #
produced by	
Medical	Exp. Date 02/13
Electronic	Target Value +/- Range Target Value +/- Range Target Value +/- Range
www.mes-	Automated Automated 0.0 0.0
global.com)	
Please note: QwikCheck beads verify that the concentration channel of the SQA-Vb is optimally	Report Save Cancel Enter the required information into the table. This information is on the QwikCheck™ beads box. Press Report to print out the settings. Press Save to apply the settings. Click Controls > Test Results to view the control records:
functioning	Controls > Test Results
problem	Export Report CARLE
impacts the	Controls Test Results
test results.	Number of Hecords 4
	Date Test Performed Control Level Lot Number Exp. Target Range Concentration MSC In Results (Mmp) Bange
NOTE:	2/19/2013 16:14 Level 1 011112001 10/13 45.0 6.3 47.6 NA Normal
Remember to	Click the Test Controls button and the screen below will be displayed:
set-up Controls	B-Sperm Management System
a new box of QwikCheck Beads.	Controls > Test Controls Fex XX Test Bull Met Data Met Data Controls Controls > Test Controls Fex XX Controls Met Data Controls Level Level 1 • Lot Number 011112001 Expiration Date 10/13 . Target Value +/- Range 45.0 +/- 6.3
	Continue >> Cancel
	Select the Control Level.
	• The pre-set Lot Number Expiration Date and Target Value + /- Range will
	be displayed automatically. Press Continue and follow the screen instructions
	• When testing is completed, the screen below will be displayed.
	• when testing is completed, the screen below will be displayed.
	Controls > Lest Results
	Number of Records 4
	Date Test Performed Control Level Lot Number Exp. Date Target Value Range +/- Concentration Results MSC In
	2/19/2013 16:14 Level 1 011112001 10/13 45:0 6:3 47:6 NA Normal
	 Click the Report or Export buttons to print a report or export the control results to an Excel file.

	Section 10: Ex	port		
Test Bull	Click Export button and	d two options	will be displayed:	
Test Data	I	Export		
Video		Export Test Resu	Its Export Controls	← BACK
Controls				
Export	Activation Export Tes	t Results fea	ture allows transferr	ing the entire archive of the test
Set-Un	results to the Excel file	(the user sho	ould specify the path	with the name of the file).
Exit	Activation Export Con Excel file (the user sho	trols feature uld specify th	allows transferring t he path with the nam	the entire control archive to the new of the file).
Set-Up	Section 11: Se	t-up		
	Set-up section contair	• Is 7 sub-meni	us:	
	Set-up			
	BSE Report Qwikclick	Defaults	Language Password	Port SQA - Yb + BACK
BSE Report setting	Click BSE Report to d	isplay the set	-up grid for the BSE	Report (shown below).
j	Set-up > BSE Re	port		
	BSE Report Phy	sical Examination	Testing Facility	_
	Body Condition Sco Body Condition Type			Show This information in Report
	Pelvic Height Pelvic Width		Address	Logo
	Pelvic Area Feet/Legs			Browse
	Eyes Vesicular Glands			Electronic signature
	Ampullae/Prostate		City	Browse
	Penis/Prepuce		Zip	
	Epididymides		State	
	Scrotum (Shape) Other]
	Scrotal Circumferen	ce 🔽 nen Examination	Phone Number	7
	Collection Method Response		Fax Number	7
	Serr Concentration (M/ml	en Characteristics	E-Mail	
	Motility (%) Progressive Motility	(%) V		
	% Normal Cells % Primery Abnorma		Web Site	7
	% Secondary Abno	rmalities		1
	WBC RBC			
	Other			
	☑ select/deselect a	II		
	Apply	Cancel	l	
	The Physical and Sen parameters recommen deselected (by clicking	ten Examina ded by the So) so that they	ition parameters dis ociety for Theriogence will not appear on t	played in the BSE report are blogy. Parameters can be the BSE report.
	Press: Apply to save t	he settings or	Cancel to not save	the entries.

OwikClick	
setting	Select: QwikClick ¹ Setting: Yes / No to begin.
	each time after running a test.
	• If No is selected the QwikClick screen will not be displayed and the automated
	morphology will be reported. (QwikClick Mcan be opened off-line to perform a morphology assessment. These results will then be displayed in the BSE report.
	Set Morphology abnormalities using the following table:
	Classification Assign Key Hep
	Primary Abnormatities Underdeveloped F2 V Use this table to:
	Double Forms F3 × Heads F4 × Memory F5 •
	Froctional Droplets F6 Trais F7 -
	Secondary Abnormalities all dividual mod. Heads F8 • Absolution mod. • The NORMAL cell is locked and cannot be changed.
	Dittal Droplets F10 Z. ASSIGN KEY (FN): Tails F11 Click on each cell in the ASSIGN KEY column.
	Other F12 Apply Cancel • Select a FUNCTION KEY (FN) from the drop down menu in each cell. • The NORMAL key is locked - it is always F1.
	The classification of morphology defects and function (F) keys for counting the different
	morphological abnormalities are manufacturer defaults. They can be re-set by the user ner Help instructions
	Defende Click this butter and the sub-many shown below will be enabled
	Defaults - Click this button and the sub-ment shown below will be opened.
Default setting	Set-up > Defaults
	Cut-off Settings Normal Ranges Concentration Std. Failure Lodes Archive Management CBALK
Cut-off settings	Cut-off Settings can be pre-set by entering (overwriting) the cut-off values in the table provided:
	Set-up > Defaults > Cut-Off Settings
	Cut-off Defaults
	Sperm Conc. [M/ml]
	Mobility [%] 40.0 Prog. Mobility [%] 30.0
	Normal Lells [%] /U.U
	Apply Cancel
	• Sperm Concentration Cut-off (M/ml): Enter the lowest acceptable Concentration (leave it empty if there is no cut-off)
	 Motility Cut-off (%): Enter the lowest acceptable Motility.
	Progressive Motility Cut-off (%): Enter the lowest acceptable Progressive
	Motility.
	• Normal Cens Cut-on (%): Enter the lowest acceptable % Normal Cens.

Normal Ranges setting	Normal Ranges -	 Set-up normal r Set-up > Defaults > 	anges for sem Normal Ranges	ien test re	sults:		
j	Normal Ranges 🗲 BACK						
		Semen Parameters		lormal Range			
		Concentration [M/ml]		100.0 - 2000.0	(
		Motility [%]		30.0 - 100.0	(
		Progressive Motility [%]	25.0 - 100.0	(
		Normal Cells [%]		70.0 - 100.0	(
		Primary Abnormalities [[%]		(=		
		Secondary Abnormaliti	ies [%]		(
	 Click on th Normal Ra BETWEEN, 	e ARROW + wh ange Settings s etc.) and the MIN	nich points to screen will ap I/MAX values,	the para opear. Er press Apr	meter to nter a syr bly .	be set-up. The nbol (<u><</u> , ≥, =,	
		Normal Range Settings					
		Field Concentration [M/ml]	Symbo BETWEEN	Min 100.0	Max 2000.0		
		Apply C	Cancel				
Conc. Standard	Concentration St Standard 1 correst	andard – Click to onds to Hemocvt	this button to	set-up Co ing chamb	oncentration oncentration on the second s	on Standard 1 or 2. o Nucleocounter.	
setting	ľ	,		5			
Failure Codes	Failure Codes – S	Set the failure cod	les that best i	ndicate the	e reason tł	ne sample failed.	
setting	Set-up > Defaults	> Failure Codes					
	# Description of Failure						
	2 Low N	lorphology lotility					
	3 Low C	oncentration					
Archive	Archive Manager	nent: Click this b	utton to view	the screer	n below. I	t is best (because	
setting	of the LARGE data	base of bull infor	mation) to set	t records t	o be archiv	ved in MONTHS	
-	(rather than years	 and as short a 	timeframe as	possible).			
	Entor						
	Enter. Archis	o records older then:	Mar	othe			
	Archive records older than: 6 Months						
				A	pply (Cancel	
Language	Click: Language	to customize the	language used	l in B-Spe	rm Report	s. Choose "Other"	
setting	from the Language	e drop-down men	u, fill-in the ta	ble as des	sired and c	lick Apply .	
		Set-up > Language	-				
		Language C BACK					
		# Bulls	# Bulls	E	English 👻		
		# Motile Sperm - Frozen # Motile Sperm - Frozen	# Motile Sperm				
		# of Dilution Steps # Prog. Motile Sperm - Fresh	# of Dilution Steps # Prog. Motile Sperm				
		# Prog. Motile Sperm - Frozen # Sperm - Fresh	# Prog. Motile Sperm # Sperm				
		# Sperm - Frozen Abaxial Implantation [%]	# Sperm Abaxial Implantation [%]				
		Abnormalities Address	Abnormalities Address				
		Age Ava. PMLF [%]	Age Avg. PMLF [%]				
		Avg. PMLF Cut-off	Avg. PMLF Cut-off				
			A DMLET (
		Apply Cancel					

Password setting	Click Password to the screen below. Cl Set-up > Password Language ← BACK	change to a r ick Apply to i	iew pa mmed	ssw ate	ord. Enter ar ly re-set and	nd conf save t	firm the new password in he new password.
	New password Confirm new password						
	Apply Cance	el					
Port setting	Port: Click this butto	on to set the o	commu	nica	ation port for	the PC	2:
		Set-up > Port					
		Password	•	BA	CK		
		Communication	port —				
		⊙ COM1	O COM	2	O COM3 C) COM4	
		O COM5	O COM	6	O COM7 C) COM8	
		О СОМЭ	О СОМ	10	O COM11 C) COM12	2
		Apply		Ca	incel		
SQA-Vb setting	SQA-Vb – Click this	button and tw	vo opti	ons	are displayed	d:	
	 Self-Test Dat Verification 	а					
	Click Self-Test Data	a and this info	ormatic	n w	vill be displaye	ed:	
		Motility Ch	annel		Conc. Cl	hannel	
		REF. I	170	mv	REF. Z	2800	mv
		AMPLITUDE	o 67	m∨	OD 1	0.0	M/ml
		ZERO LEVEL	512		OD 2	1.2	M/ml
					OD 3	2.4	M/ml
			A	gorit	hm		
		OD	0.0		CONC.	0.0	M/ml
		COUNT	33.5		MOTILITY	0.0	% M/ml
		AVERAGE WIDTH	15712.0		MOTILITY GRADE	1.0	14171111
		NUMBER SPIKE	63		MORPHOLOGY	NA	%

SQA-Vb	Print the Service Data Re	port v	which contair	ns ado	litiona	li	nformation	:
setting	Zoom	75%						
			Ser	vice D	ata Re	эр	ort	
			SQA-Vp SN		633			
		P	rint Date & Time	•	11/08/20	01 [.]	1 18:02:08	
			Self-Test Data				Service	e Data
	Ref. 1		170 [170]	m∨		1	Control Ref. 1	11
	LED Cu	r.1	8 [8]	mA		2	Min. Sp. Height	5
	Amplitu	de	67	m∨		3	Max. Sp. Width	150
	Zero Le	vel	512			4	SMI Thresh.	27
	Ref. 2		2830 [2800]	m∨		5	Average	33.52
	LED Cu	r. 2	19 [19]	mA		6	Zero Level	512
	OD 1		0.000			7	OD	0.000
	OD 2		1.215			8	MSC Amplif.	50
	OD 3		2.367			9	Min. Sp. Width	10
					- 1	10	Noise Thresh.	6
			Algorithm			11	Control Z.L.	109
	OD		0.000			12	Count	31
	Avera	je	33.52			13	OD Amplif.	97.0
	Count		31			14	Transm.	100
	Avera	ge Width	15712			15	OD Value	1.700
	Numbe	r Spike	63			16	OD Correction	110.0
	Conc.		0.0	M/ml		17	Test Noise	1
	Motility		0.0	%		18	LB OD Amp.	97
	MSC		0.0	M/ml		19	Amp. Correction	200
	Motility	Grading	1			20	Amplitude Amp.	100

The **Verification** function is for Service Personnel only.

Section 12: Exit

Click **Exit** to close the B-Sperm program. Confirm **YES** or **NO**.

Appendix I: Semen Sample Preparation

EQUIPMENT REQUIRED:

- Testing Media: QwikCheck™ Diluent for FRESH BULL semen
- Diluent Dispenser
- 10 ml Plastic Containers
- 25-micron filter
- 1 ml Pipette with tips
- SQA-Vb capillary
- Heating stage

FRESH SEMEN SAMPLES:

DILUENT/TESTING CAPILLARY PREPARATION (prior to testing the sample)

- Pre-heat testing capillaries to 37^oC / 98.6^oF for at least four minutes.
- Bring QwikCheck™ Diluent for FRESH BULL semen to room temperature (RT: 22°C-26°C or 72°F-79°F).
- Extract 2.0 ml of **QwikCheck™ Diluent for FRESH BULL semen** into a 10 ml plastic container using the dispenser (Fig. 1).
- Place 10 ml plastic container with the diluent into the heating device (see appendix section for how to operate). Pre-heat the QwikCheck™ Diluent for FRESH BULL semen at 37°C / 98.6°F for at least 4 minutes.

FRESH SAMPLE PREPARATION

- 1. Mix Fresh semen sample thoroughly.
- 2. Extract exactly 200 µl of semen using a pipette (Fig. 2).
- Add the semen from the pipette to the 2.0 ml of QwikCheck[™] Diluent for FRESH BULL semen in the plastic container (Fig 3).
- Gently but thoroughly mix the sample for 10-20 seconds (Fig 4) to insure that the diluent and semen are homogeneously mixed.
- 5. Unplug 1 ml syringe from the capillary and fill it with 1 ml of diluted semen.
- 6. Pour semen from the syringe onto the 25-micron filter.
- Open the filter cup, fill and wipe the testing capillary (preheated to 37^oC) per the instructions in the Appendix II section of this SQA-Vb User Guide.
- 8. The sample is now ready for testing.

<u>Note</u>: If testing a non-diluted sample, follow the screen instructions for sample preparation.



Figure 1: Diluent Dispenser



Figure 2: Sample Preparation



Figure 3: Add semen to the diluent



Figure 4: Mixing the sample

Appendix II: Capillary Filling Instructions

Motility Section

- Push the syringe piston in fully. Place only the thin part of the capillary into the bottom of the diluted sample -Figure 1.
- Placing two fingers below the piston head pull the piston back slowly while keeping the tip of the capillary well below the sample level and below any surface bubbles (Figure 1). Continue to aspirate the sample until it appears in the Luer adaptor.
- 3. Hold the capillary in a vertical position and visually confirm that the sample has completely filled the thin section and the cuvette section and appears in the Luer adaptor (Figure 2).
- 4. Tap on the syringe to make sure there are no air bubbles in the sample.
- 5. Quickly and thoroughly wipe both the top and bottom of the outer surface of the capillary with a tissue such as Kimwipes, etc. (Figure 3). Visually confirm that the capillary chambers are still full after wiping. If some of the sample has been depleted, a meniscus will be visible in the thin section of the capillary. If this is evident, push very slightly on the piston to re-fill the thin capillary section.
- 6. Slowly and carefully push-in the separating valve until it is level with the plastic. The capillary is now ready for testing (Figure 4)
- 7. Insert the capillary into the SQA-Vb (Figure 5)



Figure 1



Figure 2



Figure 3





Figure 5

25

NOTE: If air bubbles are still present in the capillary after

tapping on the

syringe, dip the

capillary into the

semen sample again and

aspirate a small

air bubbles into

quantity of semen to draw

the syringe.

SQA-Vb Testing Capillary

NOTE: It is important to remove all semen from the exterior of the capillary in order to prevent the SQA-Vb optical chamber from becoming clogged.

Appendix III: SQA-Vb Cleaning Instructions

When to clean:

- Daily or after every 25 tests
- If the system fails SELF-TEST

Cleaning kit components:

- Blue Dot capillaries (Fig. 1)
- Sponge-tipped drying capillaries (Fig. 2)
- Cleaning brush-wooden-handled (Fig. 4)
- Cleaning fluid

CLEANING: STEP 1

- 1. **TURN OFF** SQA-Vb and unplug it at main electrical outlet.
- 2. Select a **BLUE DOT cleaning** capillary (Fig. 1).
 - Moisten with ONE drop of cleaning fluid, shaking off excess fluid.
 - Insert into the measurement compartment fibrous material facing up, and move back and forth a few times in the directional runner.
 - Repeat with fibrous material facing down
 - Select a sponge material capillary (Fig. 2) and insert it in the same compartment in order to dry the chamber (Fig. 3)

CLEANING: STEP II

Clean the channel that measures concentration using the cleaning brush (Fig. 4):

- 1. Insert the brush (bristle-side down) fully into the upper portion of the lower chamber of the SQA-Vb in same manner as a testing capillary (Fig. 5).
- Pull the brush out of the chamber while sweeping or "dusting off" the LED (you will feel a step or shelf at the back and top of the chamber – this is the top of the LED) (Fig. 6).
- 3. Switch SQA-Vb unit **ON** and observe self-test results. The SQA-Vb should now PASS the self-test. If not, repeat cleaning procedure with the brush.











Appendix IV: Capillary Washing/Drying Instructions

(For animal applications ONLY!) Both testing capillaries and 10ml sample collection cups can be washed and re-used up to 5 times by following this EASY procedure:

Testing Capillary



Repositioning the blue valve with the jig



Removing the plunger



Reassembled capillary

Washing 10 ml sample collection cups

Step 1 After running a test:

- Use the white capillary jig to re-position the blue capillary valve
- Expel semen by pumping the plunger a couple of times
- Soak the testing capillary in tap water until ready to wash

Step 2 Set-up: Fill with 1 liter/2 quarts of solution as follows:

- Bowl #1: Tap water (marked "TAP WATER")
- Bowl #2: Distilled water (marked "DISTILLED WATER")
- Bowl #3: Isopropyl Alcohol 70% 100%

Step 3: Remove all liquid from the testing capillary:

• Pump the syringe plunger a couple of times to remove liquids.

Step 4: Capillary Washing - Follow this order:

- Bowl #1 Tap Water: Completely fill each capillary with tap water. Expel the solution into a hazardous waste container. **Repeat 2 times** then go to Bowl 2.
- Bowl #2 Distilled Water: Completely fill each capillary with distilled water. Expel the solution into a hazardous waste container. Repeat 2 times then go to Bowl 3.
- Bowl #3 Isopropyl Alcohol: Completely fill each capillary with isopropyl alcohol and expel the solution into a hazardous waste container. Repeat 2 times.
- Remove the plunger from the syringe.

Step 5: Capillary Drying Options:

- Place the capillaries:
 - On a flat surface and dry overnight.
 - o In a commercial desiccator follow manufacturer instructions.
 - In an oven on low heat for a few hours.

Step 6: Final Preparation/Inspection:

- Replace the plunger into the syringe and inspect the capillary.
- Discard capillaries with debris, cracks or broken parts.
- Make a dot on the capillary with a water proof marker after each washing cycle to indicate the # washings.

Capillary re-assembly

- Confirm (using the jig) that the blue stopper is in the correct position.
- Check syringe by aspirating air in and out twice.

Washing 10 ml sample collection cups

Refer to Step 4 and Step 5 of the Capillary Washing Procedure above - follow the same process for washing in solution bowls #1; #2 and #3. Turn upside down on absorbent paper to dry overnight or place in a commercial warming oven for a few hours.

NOTE: The focus knob should not be forced to turn beyond the initial stopping point. Turn the knob gently and when resistance is felt it is at the maximum (or minimum) position. Forcing this knob beyond the stopping point will cause extensive damage to the SQA-Vb.

Appendix V: The Visualization System

The SQA-Vb Visualization System permits the user to analyze/view semen samples using either a standard slide or the SQA-Vb testing capillary. Additionally, the visualization system is a critical "link" to B-Sperm which allows the user to view samples on a PC monitor and capture and save video images. The visualization system:

- Accommodates both an SQA-Vb testing capillary and a standard slide.
- Control knobs located below the visualization screen can be used to set the focus, brightness and contrast.
- Operate the zoom, illumination and monitor on/off functions from the SQA-Vb keypad.
- Magnification can be changed within a range: x300 through x500 by pressing ZOOM IN or OUT.

Operating Instructions

Slide Preparation:

- Use 10 µl of semen (20 micron sample depth).
- Use only a standard slide with a 22mm x 22mm cover-slip.
- Load the prepared slide into the SQA-Vb slide adaptor and insert into the visualization compartment of the SQA-Vb.

Testing Capillary Preparation:

- Fill the SQA-Vb testing capillary following the instructions in the appendix section of this guide.
- Insert the capillary into the visualization compartment of the SQA-Vb.

Testing Process:

- 1. The video display will automatically illuminate when the SQA-Vb is turned on.
- 2. Use the monitor **ON/OFF** key on the keypad to independently operate the video display.
- 3. To ensure that the visualization system is working properly prior to use:
 - a. Press the **HIGH ILLUMINATION** key multiple times to ensure a maximum level setting.
 - b. Turn **BRIGHTNESS**, **CONTRAST** and **COLOR** buttons all the way counterclockwise.
 - c. Turn FOCUS knob fully clockwise.
- 4. Use **ZOOM IN** for maximum magnification (x500) and **ZOOM OUT** for minimum magnification (x300).
- 5. Insert the semen sample into the visualization chamber.
- 6. Turn the **BRIGHTNESS** knob clockwise until the video screen just begins to lighten-up.
- 7. Turn the **FOCUS** knob counter-clockwise until the image is in focus.
- 8. Adjust **CONTRAST, COLOR, BRIGHTNESS, FOCUS** and object **ILLUMINATION** controls for optimal image quality.

Appendix VI: Sample Heating Options (Block Heater and Stackable System) For Use with A-Tech Sperm Quality Analyzers

BLOCK HEATER

Safety:

- Do not pre-fill the SQA testing capillaries with semen prior to heating.
- Do not touch the heating elements to check the temperature.
- If transported or stored in humid conditions, dry the unit before connecting it to power.
- Plug into a grounded electrical outlet that delivers the appropriate voltage indicated on the rear panel of the heater.
- If liquid is spilled inside the unit, disconnect the power supply, take out the metal heating blocks and wipe the spilled material with a damp cloth. Do not use chemical cleaning agents.
- Place the system on a level surface that is free from flammable materials, insuring that all ventilation slots on the base of the system are clear of obstructions.



Heater Operation:

- Place the metal heating racks for capillaries and sample containers into the heating system.
- Plug the electrical cord into the socket at the rear of the unit.
- Turn the unit on by pressing the Main Switch on the front panel The display will illuminate showing the current temperature of the block.
- Press the 'S' button and adjust the temperature to 40°C using the Control knob. This will ensure a 37°C temperature for the testing capillary and semen samples.
- Press 'S' to confirm the temperature setting or press 'F' to exit without changing the value.
- The block heater will now begin to heat the racks.
- A light will indicate that the system is heating and will begin to flash when the set temperature is approached.
- The temperature setting will be stored in the memory.

Timer Operation:

- Press 'F' button and turn the Control knob until 'CLOC' is seen.
- Press 'S' button and turn the Control knob until 'On' is seen.
- Press 'S' button and turn the Control knob to select the required time (Example: 4 minutes (004) for Bull semen).
- Press 'S' button: The time and the Temperature will be displayed intermittently.
- A beep will sound and 'End' will be displayed when the time has expired; press 'S' to stop beeping.

Heating the SQA Testing Capillaries and Samples:

- Place **<u>empty</u>** testing capillaries in the heating rack (as shown in the picture above).
- Place empty 10-ml plastic containers in the appropriate heating rack.
- Wait 5-7 minutes for the heating unit to pre-heat the containers and testing capillaries.
- Distribute the diluent and add the semen sample into 10-ml container following User Guide instructions.
- Close the containers during pre-heating.
- Remove testing capillaries as needed for testing. Fill with semen per user guide instructions.

Note: If the processor detects an error in heating the fault indicator will illuminate, the buzzer will beep and the display will flash. To reset this fault, switch the unit off and on. If the fault re-occurs contact the service personnel at your local distributor.

Appendix VII: Glossary of Terms

	SQA-Vb Terms	Definition
	SN	Serial Number of the SQA-Vb
	DATE/TIME	The date and time the test was performed
Sample/Test Data	SAMPLE #	The number assigned to the semen sample
	BULL ID	The identifying number of the bull being tested
	BREED	Identification of the bull breed
Test Results	CONCENTRATION	Sperm Concentration expressed in millions/ml
	MOTILITY	% of Motile Sperm
	PROGRESSIVE MOT.	% of Progressively Motile Sperm
	NORMAL CELLS	% of Morphologically Normal Sperm

Appendix VIII: SQA-Vb System Specifications

Dimensions:	40 x 30 x 15 cm
Weight:	4 kg
AC power supply:	100 to 250 VAC, 50/60 Hz, 10 VA

Measurement Compartment

- Sources of radiant energy two 880 nm LEDs for motility and spectrophotometry channels
- Detector system 2 photo detectors Motility and Optical Density

Visualization Compartment

- Green LED illumination system
- CCD, 350 TV lines
- Objective lens: Standard, x20
- Signal Output: PAL standard
- Zoom system for smooth magnification transition from x300 to x500
- Focus regulator

Display(s)

- Operational backlight LCD (16 lines x 40 characters)
- Video backlight LCD (8 x 10 cm)

Printer

- Built-in, Dot Matrix
- Non-thermostatic narrow paper with 20 characters per line (Citizen)
- Ribbon cassette (Citizen)

Keypad

- **Operational keys:** ON/OFF, TEST, PRINT, SERVICE, ARCHIVE, DELETE, ENTER, four cursor buttons, ESC, numeric buttons (0-9)
- Video control keys: ZOOM IN/OUT, ILLUMINATION HIGH/LOW, and MONITOR ON/OFF

Front Panel

- Built-in printer
- Visualization compartment
- LCD video display and controls
- Focus knob
- LCD operational display
- Measurement compartment
- Multi-button keypad

Rear/Side Panel

- Power connector with fuse-holder (fuse 250V, 1A)
- Video connector
- RS232 cable outlet
- I-Button port (side panel)

Specimen Testing Supplies

- **Measurement capillary:** Disposable, multi use plastic, positive displacement testing capillary (purchase from manufacturer).
- Standard lab slide: 76 X 25.6 mm with 22 x 22 mm cover-slip.
- I-Button: Required to run tests (purchase from manufacturer)

Archive Capacity

• 500 test records in each archive (Test and Controls)

Operating System

- Control: Keypad
- Analysis Time: Normal Test 50 seconds
- **Software:** Resides on flash memory and drives all man-machine interface functions, runs algorithms for test measurements and operational screens. System can be upgraded from a PC CD-ROM.
- **Sample Testing Temperature:** Calibrated for room temperature only. Motility results will be impacted by over-heating the specimen.
- Motility channel input signal: Analog, up to 5V.
- Spectrophotometer channel input signal: Modulated (1 kHz) analog, up to 5V.

Quality Control

- Internal: Electronic Self-Test and Auto-Calibration.
- **External:** Latex beads control material (QwikCheck[™] beads)

PC Compatibility

Minimum requirements for B-Sperm[™] software

- PC: 1 GHz processor, Pentium 3
- **RAM:** 256 MB
- AGP-video display card with at least 16 MB of RAM memory
- Video color: At least 16 bit (65,535)
- CD ROM drive
- 40 GB free hard disk space for image capturing/storage
- Video resolution: Minimum 640 x 480

Operating system compatibility

- Windows XP or VISTA or WIN7
- Ports: One serial; two USB ports
- Monitor: 15" color

Additional Software (supplied with system)

• **B-Sperm software:** Real time visualization interface between PC and SQA-Vb visualization system, data transfer, video/picture capture and archive.

Operational Temperature and Humidity

- System is operational at 15-38°C.
- NOTE: SQA-Vb operates in a wide range of ambient temperatures however the system is calibrated to measure semen samples pre-heated to 37^oC / 98.6^oF prior to testing (refer to the SQA-Vb User Guide Appendix I: Semen Sample Preparation).
- *NOTE*: Variations in ambient temperature may impact the accuracy of test results because of the effect of temperature on semen.
- System is fully operational at up to 80% humidity and 31°C.

Maintenance Schedule

• Cleaning daily or after every 25 tests (refer to the SQA-Vb User Guide Appendix IV: SQA-Vb Cleaning Instructions).

Manufacturer Recommendations

- Operate the SQA-Vb away from devices that may cause electronic noise (cell phones) or other devices causing vibrations such as centrifuges.
- Turn system **OFF** at the rear-panel when not in use for extended period of time.
- Variations in ambient temperature can affect semen samples. The SQA-Vb is calibrated to conduct tests at 37^oC / 98.6^oF.
- Semen is considered a biologically hazardous material and is subject to individual laboratory protocols for handling such materials.

Factory Default Settings

Date format: **DD/MM/YY** Date/Time: Manufacturer's local date/time Sample Type: **FRESH** CONC STANDARD: **1 (Neubauer)** Operational Default Settings:

- Automatically Send to PC: YES
- Data Entry in B-Sperm: YES
- QwikClick Counter: NO

Appendix IX: Product Performance Data

Abbreviations:

Sperm Concentration
Motile Sperm Concentration
Progressively Motile Sperm Concentration
Micron
Second
Million per milliliter

Performance Data Summary:

The performance of the SQA-Vb system for bull semen analysis is summarized in the text, tables and graphs below. Sperm concentration measurements are expressed as x10⁶ sperm cells per milliliter (M/ml). Motility and Morphology values are expressed as a percent (%). Unless otherwise noted all testing was performed using fresh and frozen bull semen samples.

Calibration:

Each SQA-Vb device is biologically calibrated against two reference systems at Medical Electronic System's laboratory using bull semen.

Dynamic Range:

Sample Type	Test Mode	Conc. M/ml	Motility %	Morphology %	MSC M/ml	PMSC M/ml	Velocity, mic./sec.
Fresh	Fresh	0-2000	0-95	0-100	0-1900	0-1800	0-130
Frozen	Frozen	0-200	0-95	-	0-200	0-200	0-80

Precision and accuracy established against a known target (Latex beads)

Background: The precision and accuracy of the SQA-V was compared to a known target value using commercially available latex beads of two concentrations. Latex beads are run on the SQA-V in the same manner as semen samples.

Limitations of method:

Latex beads cannot:

- Measure sperm motility or morphology
- Correct for inaccurate chamber depths or technician errors

Method comparison:

A total of 320 latex bead samples were tested on ten SQA-V systems (32 samples per SQA-V). Precision of the SQA-V was estimated (Table 1). SQA-V concentration readings were compared to the established target values +/- the acceptable range for the latex beads (Fig. 1 & 2).

Accu-beads® published acceptable ranges

- (Hemacytometer):
- Vial #1: 47 +/- 7.0 M/ml
- Vial #2: 24 +/- 3.4 M/ml

Table 1: Precision

SQA-V	Accu-beads®	CV %
Intra-device	High 47±7.0 M/ml	≤ 0.01
Variability	Low 24 ± 3.4 M/ml	≤ 0.01
Inter-device	High 47±7.0 M/ml	≤ 2.00
Variability	Low 24 ± 3.4 M/ml	≤ 2.50

Fig. 1- 2 Accuracy Low/High Level Controls:





Conclusions:

The CONTROL mode software of the SQA-Vb (bull) device is exactly the same as the SQA-V (human) system. Both systems also have the same hardware platform. Therefore, the accuracy and precision results obtained on the CONTROL mode of the SQA-V will be the same as that of the CONTROL mode of the SQA-Vb.

Sensitivity, specificity, precision, accuracy and correlation to manual method established in the MES laboratory and field clinical trials using bull semen samples

FRESH SEMEN Performance Claims (Table 2-5)

Sensitivity

- Concentration: 90%
- Motility: 85%
- Morphology 80%

Specificity

- Concentration: 90%
- Motility: 80%
- Morphology 90%

Precision (CVs)

- Conc.: 3%
- Motility: 5%
- MSC: 7%
- Velocity: 10%

Accuracy (regression coefficients of "live/dead" trend line) Claims:

- Motility: 0.85
- MSC: 0.9
- PMSC: 0.9

Correlation to Manual Method Claims:

- Concentration: 0.9
- Motility: 0.8
- Morphology 0.7
- MSC: 0.9
- PMSC: 0.8
- Velocity: 0.75

FROZEN SEMEN Performance Claims (Table 3, 5)

Precision (CVs)

- MSC: 7%
- Velocity: 4%

Correlation to Manual Method

- MSC: 0.8
- PMSC: 0.7
- Velocity: 0.85

Table 2: Sensitivity/Specificity

SQA-Vb vs. Microscope	Sensitivity	Specificity	
FRE	SH SEMEN		
Sperm Concentration M/ml	100.0%	98.4%	
Motility, %	96.3%	85.0%	
Morphology	85.7%	91.7%	

Table 3: Precision SQA-Vb intra-device Variability (CV, %)

	Sample Type			
Parameter	Fresh	Frozen		
Sperm Concentration M/ml	2.4	-		
Motility, %	4.1	-		
Morphology %	5.0	-		
MSC, M/ml	6.1	6.4		
Velocity, microns/sec.	9.9	3.4		

Table 4: Accuracy: Regression coefficients from "live/dead" experiments

Parameter	SQA-Vb	Manual
Motility, %	0.9135	0.8325
MSC, M/ml	0.9348	0.7547
PMSC, M/ml	0.9340	-

Notes:

- Sensitivity and specificity claims are lower than actual values noted (Table 2).
- Precision CV claims are higher (lower precision) than actual values noted (Table 3).
- Accuracy regression coefficient claims are less than actual values noted (Table 4).
- Correlation to Manual Method claims are less than actual correlations noted (Table 5).

Method comparison:

SQA-Vb was compared to the microscope based on WHO '99 guidelines. The SQA-Vb automated readings of the sperm concentration, motility, MSC, PMSC and velocity were compared to microscopic results. A Makler chamber was used according to manufacturer's instructions for the manual sperm concentration measurements. A standard slide and B-Sperm software were used to assess manual motility, progressive motility and velocity measurements. Manual MSC and parameters were calculated from PMSC experimental results. The protocols were based on WHO '99 and MES guidelines. The alpha-site clinical trials were conducted at the Sion farm. A total of 104 fresh and 138 frozen semen samples were analyzed.

Analytical Specificity:

- To achieve analytical specificity a specific wave length of light which is maximally absorbed by sperm cells and minimally absorbed by other cells and seminal plasma is used.
- Low noise and high electronic resolution hardware components and compensation circuits ensure that analytical specificity is optimized.

Limitations of method:

Samples were assessed in duplicate on the automated SQA-Vb system and manually using a microscope. Statistical counting errors and intra-operator variability (subjectivity) may have affected the results of the study.

Accuracy assessment: "live/dead" sperm plots.

The SQA-Vb accuracy was assessed by "live/dead" bull sperm experiments. Fresh bull semen was distributed into two aliquots. The first aliquot was intact ("live") and the second one was treated with the liquid nitrogen ("dead"). Then different "live-to-dead" proportions were created providing constant Sperm а Concentration but varying MSC, PMSC and Motility. The samples were tested using the SQA-Vb device and under the microscope and the results plotted. The linear trend lines were established for motility, MSC and PMSC variables vs. "Live/Dead" sperm ratio.

Table 5: Correlation to manual method

Deremetere	Correlation coefficients				
Parameters	Fresh semen	Frozen semen			
Sperm Concentration, M/ml	0.93	-			
Motility, %	0.81	-			
Morphology %	0.71	-			
MSC, M/ml	0.94	0.84			
PMSC, M/ml	0.86	0.74			
Velocity, mic./sec.	0.81	0.91			

Fig. 3: Method comparison: Regression plot of SQA-Vb Sperm Concentration in fresh bull semen vs. manual results



Fig. 4: Method comparison: Regression plot of SQA-Vb Motility in fresh bull semen vs. manual results



Fig. 5: Method comparison: Regression plot of SQA-Vb MSC in fresh bull semen vs. manual results



Performance parameters:

- Sensitivity and specificity were calculated using ROC analysis formulas. The cutoffs normally used for the sperm concentration and motility of the fresh semen samples were used for calculation of sensitivity and specificity. As for the frozen samples there are no cutoffs, sensitivity and specificity were not calculated for this type of samples (Table 2).
- Precision of the SQA-Vb device was estimated by calculation of the intra-device coefficients of variation (CV) of the duplicate measurements (Table 3). CV is calculated according to the formula:
 CV = SD / MEAN x 100
- The accuracy of the SQA-Vb device was characterized by regression coefficients of the trendline obtained in the "live/dead" experiment (Table 4).
- Correlation to manual method was established by calculation of correlation coefficients (Table 5, Fig. 3-5).

Conclusions:

- High levels of sensitivity, specificity and correlation of the SQA-Vb device results to the manual method were found. Therefore the instrument can be used in the field for the semen quality assessment, dose preparation and in the frozen semen QC.
- SQA-Vb provides the precise and accurate results with low coefficients of variation (<6%) and high regression coefficients of the "live/dead" trend lines (>0.95).

	Appendix X: SQA-Vb Service Menu			
	Press Service button of the SQA-Vb system to open SERVICE MENU and access one of 5 functions:			
	SERVICE DATA: Select this option to view three service screens:			
	• Service Data screen (Communication screen): Establishes communication between the SQA-Vb and the PC (B-Sperm) in order to:			
	 Set-up the CONTROLS default settings Report SQA-Vb service information for technical troubleshooting Self-Test Data screen (1): Provides system information after Self-Test (Self-Test Data) and after regular testing (Internal Data). Self-Test Data screen (2): This screen displays internal parameters and tests remaining. 			
	 SERVICE PERSONNEL: For technical services personnel only and requires a password to access. 			
	PRINT SELF-TEST DATA:			
	Select to print the SQA-Vb Self-Test Data.			
	DATE SETTING: Select this option to set the date:			
	DATE SETTING			
	LOCAL TIME: 08:15:45 HH:MM:SS 24 h DATE FORMAT: MM/DD/YY D/MM/YY DATE SETTING: 04/01/07			
NOTE: If the I-Button is not properly inserted a message: I-BUTTON NOT PROPERLY ACTIVATED will be displayed. Remove the button, press ESC and try again.	LOCAL TIME: enter local time			
	 DATE FORMAT: Select the format DD/MM/YY or MM/DD/YY using the right/left arrows on the keyboard. Press ENTER to confirm 			
	DATE SETTING: enter current date.			
	ADD TESTS TO COUNTER: I-button tests			
	Select SERVICE MENU > ADD TESTS TO COUNTER and press ENTER			
	The SQA-Vb screen will instruct the user to:			
	TO ADD MORE TESTS			
	HOLD NEW I-BUTTON AGAINST PORT			
	AND PRESS ENTER			
	 Make sure the I-Button touches the internal surface and edges of the port. Press the I- button firmly in the port, moving it left and right to make sure it also touches the edges of the port. 			
	• The # TESTS ADDED and the cumulative # TESTS NOW REMAINING will be displayed on the screen of the SQA-Vb.			
	• The screen will warn the user if an EMPTY I-BUTTON (was) INSERTED			

Appendix XI: Troub	leshooting SQA-Vb Warning Screens
Stabilization Failed:	STABILIZATION FAILED TURN OFF MAIN SWITCH ON REAR PANEL REACTIVATE UNIT
	IF PROBLEM PERSISTS, CALL FOR TECHNICAL SUPPORT
• Ensure there is no te	sting capillary in the measurement compartment.
Remove the SQA-Vb	from sources of electronic noise (centrifuge, cell phones, etc.).
Clean measurement	compartment (refer to Appendix).
• Reboot the SQA-Vb v	without a testing capillary in the chamber:
o Turn system OFF	then back ON at the main switch on the rear panel.
o Press the front pa	anel ON/OFF key to begin Auto-Calibration/Self-Test
o Call for technical	support if failure recurs.
Failed Self-Test:	FAILED SELF-TEST TURN OFF MAIN SWITCH ON REAR PANEL CLEAN OPTICAL CHAMBER REACTIVATE UNIT
	IF PROBLEM PERSISTS,
Ensure there is no te	esting capillary in the measurement compartment.
Remove the SQA-Vt	o from sources of electronic noise (centrifuge, etc.).
Clean measurement	compartment (refer to Appendix).
Reboot the SQA-Vb	without a testing capillary in the chamber:
o From the rear par	nel switch, turn the system OFF then back ON.
 Press the front pa Test. 	anel ON/OFF key to begin Auto-Calibration, Stabilization and Self-
Call for technical sup by printing a copy of	oport if this message is displayed again. Prepare for technical support the SQA-Vb Self-Test/Service Data:
o Press SERVICE I	key. The SERVICE MENU will be displayed.
o Select PRINT SE	RVICE DATA option.
o Press ENTER to	print a copy of the self-test/service data.
lectronic Noise:	ELECTRONIC NOISE TURN OFF MAIN SWITCH ON REAR PANEL REACTIVATE UNIT
	IF PROBLEM PERSISTS, CALL FOR TECHNICAL SUPPORT
• Ensure there is no te	esting capillary in the measurement compartment.
Remove SQA-Vb fro	om sources of electronic noise (centrifuge, etc.).
Clean measurement	compartment (refer to Appendix).
o After cleaning, ac	tivate MAIN MENU > TEST NEW SAMPLE and re-run the test.
• If this message is di	splayed again, report the SQA-Vb

• Turn the system OFF then back ON at the main switch on the rear panel.
 Press the front panel ON/OFF key to begin Auto-Calibration and Stabilization.
o Re-run the test.
 Call technical support if this message is displayed again. Prepare for technical support by printing a copy of the Self-Test/Service parameters:
 Press SERVICE key. The SERVICE MENU will be displayed.
o Select PRINT SERVICE DATA option.
• Press ENTER and the service parameters will be printed.
Remove Capillary:
REMOVE CAPILLARY FOLLOW ON-SCREEN INSTRUCTIONS
 If the testing capillary has been left in the measurement chamber after testing a sample the message above will be displayed.
 Remove the testing capillary before running a new test.