

MEDICAL ELECTRONIC SYSTEMS  
WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	120515	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	12-May-16	LOT# 120515001			LOT# 120515002			LOT# 120515003	
For use with:		Conc.			Conc.			Conc	MSC
SQA-V Automated		44	+/-	6.2	23	+/-	4.6	0.0	0.0
SQA-V Visualization		44	+/-	11.0	23	+/-	5.8	0.0	0.0
Hemocytometer		45	+/-	11.3	23	+/-	5.8	0.0	NA
Makler		45	+/-	11.3	23	+/-	5.8	0.0	NA
Fixed Coverslip		44	+/-	11.0	23	+/-	5.8	0.0	NA