

MEDICAL ELECTRONIC SYSTEMS

WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form											
BATCH #		010725		Level 1			Level 2			Negative Control	
EXP DATE:		July-26		LOT# 010725001			LOT# 010725002			LOT#	010725003
For use with:				Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated				47	+/-	6.6	25	+/-	5.0	0.0	0.0
Manual	SQA-V Visualization			49	+/-	12.3	26	+/-	6.5	0.0	NA
	Hemocytometer			49	+/-	12.3	25	+/-	6.3	0.0	NA
	MAKLER			49	+/-	12.3	25	+/-	6.3	0.0	NA
	Fixed Coverslip			48	+/-	12.0	25	+/-	6.3	0.0	NA