

**QwikCheck™ Low Range Beads Batch Release Form**

BATCH #	211122	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	21-Aug-23	LOT#	211122001	LOT#	211122002	LOT#	211122003		
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated		15	+/-	6.0	4	+/-	2.6	0.0	0.0
SQA-V Visualization		16	+/-	6.4	4	+/-	2.6	0.0	NA
Fixed Coverslip		16	+/-	6.4	3	+/-	2.0	0.0	NA
Hemocytometer		17	+/-	6.8	NA	+/-	NA	0.0	NA

Victoria Goldenberg, Clinical research assistant



Signature

27<sup>th</sup> November 2022

Date

Sapir Zrihan OBO Frida Cohen, Clinical Manager



Signature

27<sup>th</sup> November 2022

Date

Marcia Deutsch, CEO



Signature

27<sup>th</sup> NOV 2022

Date