



# Certificate of Analysis


Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	219634_0626
Storage Temp	Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	>= 1200	1252

Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition
-----------------	--

Expiration Date: JUN/2027

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature:   
Clinical Director

Date: 01-JUN-2026