

MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM

| QwikCheck™ Beads Batch Release Form | | | | | | | | | |
|-------------------------------------|--------------|---------|-----------|------|---------|-----------|-----|------------------|-----------|
| BATCH # | 010924 | Level 1 | | | Level 2 | | | NEGATIVE CONTROL | |
| EXP DATE: | September-25 | LOT# | 010924001 | | LOT# | 010924002 | | LOT# | 010924003 |
| For use with: | | Conc. | | | Conc. | | | Conc | MSC |
| SQA / QwikCheck gold Automated | | 48 | +/- | 6.7 | 25 | +/- | 5.0 | 0.0 | 0.0 |
| SQA-V Visualization | | 49 | +/- | 12.3 | 25 | +/- | 6.3 | 0.0 | NA |
| Hemocytometer | | 48 | +/- | 12.0 | 25 | +/- | 6.3 | 0.0 | NA |
| MAKLER | | 48 | +/- | 12.0 | 25 | +/- | 6.3 | 0.0 | NA |
| Fixed Coverslip | | 49 | +/- | 12.3 | 26 | +/- | 6.5 | 0.0 | NA |

D.K

Laboratory Technician

S.Z

Production Laboratory Manager

M.D

CEO

D.K

Signature

SZ

Signature

MD

Signature

9-Sep-24

Date

9-Sep-24

Date

9-Sep-24

Date