



## Certificate of Analysis

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	220809_0323
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	$\geq 1200$	1417.5
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: MAR 2024

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z  
Acting Lab Manager

Date: 5<sup>th</sup> Mar 2023