

MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM

QwikCheck™ Beads Batch Release Form									
BATCH #	110424	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	April-25	LOT#	110424001		LOT#	110424002		LOT#	110424003
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck Gold Automated		48	+/-	6.7	25	+/-	5.0	0.0	0.0
SQA-V Visualization		46	+/-	11.5	24	+/-	6.0	0.0	NA
Hemocytometer		50	+/-	12.5	24	+/-	6.0	0.0	NA
MAKLER		50	+/-	12.5	24	+/-	6.0	0.0	NA
Fixed Coverslip		48	+/-	12.0	25	+/-	6.3	0.0	NA

S.Z, Production Lab Lead Teach

S.Z

2-May-24

Signature

Date

M.D, CEO

MØ

2-May-24

Signature

Date