



## Certificate of Analysis

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	221113_0823
Storage Temp	Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	$\geq 1200$	1300.05
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: August 2024

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z  
Production Lab Lead Tech

Date: 1<sup>st</sup> August 2023