



Certificate of Analysis

Product Name	QwikCheck Liquefaction
Product Number	A-CA-01698-00
Product Lot	4A018929_0624
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	≥ 1200	1350

Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition
-----------------	--

Expiration Date: June 2025

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z
Production Lab Lead Tech

Date: 6th June 2024