

**MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM**

QwikCheck™ Beads Batch Release Form									
BATCH #	060325	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	March-26	LOT#	060325001		LOT#	060325002		LOT#	060325003
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck Gold Automated		48	+/-	6.7	26	+/-	5.2	0.0	0.0
SQA-V Visualization		47	+/-	11.8	25	+/-	6.3	0.0	NA
Hemocytometer		49	+/-	12.3	25	+/-	6.3	0.0	NA
MAKLER		49	+/-	12.3	25	+/-	6.3	0.0	NA
Fixed Coverslip		48	+/-	12.0	26	+/-	6.5	0.0	NA

S.L
Laboratory Technician

S.L
Signature

18-Mar-25
Date

S.Z
Production Lab Manager

S.Z
Signature

18-Mar-25
Date