



# Certificate of Analysis

Product Name                      **QwikCheck Liquefaction**  
Product Number                 0900  
Product Lot                        219634\_0426  
Storage Temp                      Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	>= 1200	1252

Unit Definition                    That amount of enzyme  
   causing a decrease in absorbance  
   at 237nm of 0.0075 per minute  
   at 25°C Unit Definition

Expiration Date:                 APR/2027

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: \_\_\_\_\_  
   *MJY*  
   Clinical Director

Date: \_\_\_\_\_  
   01-APR-2026