



Certificate of Analysis

Product Name	QwikCheck Liquefaction
Product Number	0900
Product Lot	209152_0924
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	≥ 1200	1280

Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition
-----------------	--

Expiration Date: September 2025

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z
 Production Laboratory Manager

Date: 3 September 2024