



## Certificate of Analysis

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	220809_0123
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	$\geq 1200$	1332.5
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: JAN 2024

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: SZ  
Acting Lab Manager

Date: 15<sup>th</sup> JAN 2023