



Medical Electronic Systems Ltd.  
A-Tech (a division of M.E.S.)  
20 Alon Hatavor Street, Caesarea  
Industrial Park 38900, Israel  
E-mail: [mes@mes-ltd.com](mailto:mes@mes-ltd.com)



## **Certificate of Analysis**

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	220809_1222
Storage Temp	Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	> = 1200	1332.5
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date      DEC 2023

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: Sapir Zrihan OBO Frida Cohen,  
Clinical Manager

Date: 11<sup>th</sup> December 2022

MES-0186-02