

**MEDICAL ELECTRONIC SYSTEMS  
WWW.MES-GLOBAL.COM**

QwikCheck™ Low Range Beads Batch Release Form									
BATCH #	150425	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	Jan-26	LOT#	150425001		LOT#	150425002		LOT#	150425003
For use with:		Conc.			Conc.			Conc.	MSC
SQA / QwikCheck gold Automated	15	+/-	6.0	3	+/-	2.0	0.0	0.0	
SQA-V Visualization	15	+/-	6.0	3	+/-	2.0	0.0	NA	
Fixed Coverslip	16	+/-	6.4	3	+/-	2.0	0.0	NA	
Hemocytometer	15	+/-	6.0	NA	+/-	NA	0.0	NA	

\_\_\_\_\_  
I.O, Laboratory Technician

I.O  
\_\_\_\_\_  
Signature

\_\_\_\_\_  
20-Apr-25

\_\_\_\_\_  
S.Z, Production Laboratory Manager

SZ  
\_\_\_\_\_  
Signature

\_\_\_\_\_  
20-Apr-25