

MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM

| QwikCheck™ Low Range Beads Batch Release Form | | | | | | | | | |
|---|--------|---------|-----------|-----|---------|-----------|-----|------------------|-----------|
| BATCH # | 050824 | Level 1 | | | Level 2 | | | NEGATIVE CONTROL | |
| EXP DATE: | May-25 | LOT# | 050824001 | | LOT# | 050824002 | | LOT# | 050824003 |
| For use with: | | Conc. | | | Conc. | | | Conc | MSC |
| SQA / QwikCheck gold Automated | | 16 | +/- | 6.4 | 3 | +/- | 2.0 | 0.0 | 0.0 |
| SQA-V Visualization | | 15 | +/- | 6.0 | 4 | +/- | 2.6 | 0.0 | NA |
| Fixed Coverslip | | 15 | +/- | 6.0 | 3 | +/- | 2.0 | 0.0 | NA |
| Hemocytometer | | 15 | +/- | 6.0 | NA | +/- | NA | 0.0 | NA |

Laboratory Technician

S.Z, Production Lab Lead Tech

M.D, CEO

D.C

Signature

SZ

Signature

MCD

Signature

11-Aug-24

11-Aug-24
Date

13-Aug-24
Date