



Certificate of Analysis

Product Name	QwikCheck Liquefaction
Product Number	A-CA-01698-00 / EU-A-CA-02199-00
Product Lot	214001_0725
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	≥ 1200	1247
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: July 2026

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z.
Production Laboratory Manager

Date: 16-July-25