



# Certificate of Analysis

Product Name **QwikCheck Liquefaction**  
Product Number 0900  
Product Lot 4A018929\_0824  
Storage Temp Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	$\geq 1200$	1350

Unit Definition That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition

Expiration Date: August 2025

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z  
Production Lab Lead Tech

Date: 15<sup>th</sup> August 2024