



## Certificate of Analysis

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	4A018929_0724
Storage Temp	Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	$\geq 1200$	1350

Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition
-----------------	--

Expiration Date: July 2025

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: SZ  
Production Lab Lead Tech

Date: 1<sup>st</sup> July 2024