



Certificate of Analysis

Product Name **QwikCheck Liquefaction**
Product Number 0900
Product Lot 214001_0125
Storage Temp Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	\geq 1200	1247

Unit Definition That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition

Expiration Date: January 2026

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: S.Z
Production Laboratory Manager

Date: 2 January 2025