



Certificate of Analysis


Product Name	QwikCheck Liquefaction
Product Number	A-CA-01698-00 / EU-A-CA-02199-00
Product Lot	219634_0526
Storage Temp	Room temperature

Test	SPECIFICATION	RESULTS
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	≥ 1200	1252

Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition
-----------------	--

Expiration Date: MAY/2027

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: 
Clinical Director

Date: 03-MAY-2026