

## Your abstract submission has been received

Print this page

You have submitted the following abstract to the ASRM 2024 Scientific Congress & Expo. Receipt of this notice does not guarantee that your submission was complete or free of errors.

---

### PRECISION AND ACCURACY OF THE SQA-IOW SPERM QUALITY ANALYZER IN THE HANDS OF CLIA WAIVED RN OPERATORS

---

**Natan Bar-Chama, M.D.**, Reproductive Medicine Associates of New York Department of Urology The Mount Sinai Hospital, New York, NY and Lev Rabinovitch, PhD, Medical Electronic Systems, Encino, CA

**Title:**

PRECISION AND ACCURACY OF THE SQA-IOW SPERM QUALITY ANALYZER IN THE HANDS OF CLIA WAIVED RN OPERATORS

**Submitter's E-mail Address:**

nbarchama@gmail.com

**Preferred Presentation Type:**

Oral or Poster

**Study Type:**

Case Control Study

**Category - Subcategory(ies):**

**Technology & Innovation:** Reproductive Urology

**Funding:**

This study was funded by Medical Electronic Systems (MES) Encino, California

**\* Submission of an abstract for consideration for presentation implies that the presenting author & associated co-authors have legal and ethical rights to submit and present this work. Plagiarism and submitting work that an author has no rights to, will result in an investigation and penalty.**

**\* I verify that I am in compliance with HIPAA standards to protect the privacy of the patients discussed in my presentation(s). I either have received written authorization from the patient, have removed any identifiable images or patient records from my presentation, or my presentation does not pertain to patient treatment.**

**Permissions - Prior Publication or Presentation**

This abstract contains original work, not published or presented previously at a meeting of another national or international scientific organization prior to this meeting and has not been submitted for publication at the time of this submission.

**ACCME Disclosure**

Nothing to disclose. No off-label or otherwise non-approved product use.

**Did this abstract require approval by a local Institutional Review Board (IRB) or equivalent?**

This abstract has been approved by a local Institutional Review Board (IRB) or equivalent.

**Not applying for an award****Abstract Text:**

**OBJECTIVE:** To demonstrate the precision and accuracy of the point-of-care SQA-iOw semen analyzer in the hands of non-professional waived Registered Nurse (RN) users.

**MATERIALS AND METHODS:** The new SQA-iOw is a small-footprint, point-of-care automated medical device that performs a complete semen analysis in 75 seconds. It is designed for low-

complexity Clinical Laboratory Improvement Amendments (CLIA) waived clinician offices that lack certified andrology labs or professional technicians. This study demonstrated the accuracy and precision of the Sperm Quality Analyzer (SQA-iOw) in the hands of CLIA-waived RN operators compared to the results generated by the FDA-cleared CLIA requiring SQA-V laboratory analyzer run by professional technicians.

An SQA-iOw accuracy trial was conducted at three study sites encompassing 398 donor semen samples run in singleton; once on the SQA-iOw, by nine waived RN operators, and then again on the FDA-cleared SQA-V, by professional operators. Correlations between these two groups were established for sperm concentration, motile sperm concentration (MSC), progressively motile sperm concentration (PMSC), motility, progressive motility, and sperm morphology using Passing-Bablok regression analysis.

SQA-iOw precision trial was conducted at a single site to compare the results of three CLIA-waived SQA-iOw RN operators testing seven semen samples in duplicate, spanning low, middle, and high levels of four semen parameters on three SQA-iOw devices over 0, 20, 40, and 60-minute time intervals. This resulted in 168 replicate results for four semen parameters. Statistical evaluation included mean, standard deviation (SD), and coefficient of variation (CV, %).

**RESULTS:** SQA-iOw accuracy of  $\geq 96\%$  was demonstrated for the 398 semen samples between the SQA-iOw CLIA waived RN user's and the SQA-V professional technician's results for sperm concentration, motile sperm concentration, progressively motile sperm concentration, sperm motility, progressive motility, and sperm morphology (Table 1). SQA-iOw precision was demonstrated by a  $<7.5\%$  coefficient of variation CV for sperm concentration, MSC, PMSC, and morphology measured across the 168 replicates.

**Table 1. SQA-iOw Versus SQA-V n=398**

Sperm Parameters	Intercept	Slope	Correlation
Concentration M/ml	0.05	0.98	0.99
Motility %	1.94	0.94	0.96
Progressive Motility %	-0.26	0.94	0.98
Morphology %	-1.00	1.00	0.97
MSC M/ml	0.27	0.95	0.99
PMSC M/ml	-0.17	0.92	0.99

**CONCLUSIONS:** The point-of-care SQA-iOw demonstrated a high level of accuracy and precision in the hands of CLIA-waived RN users.

**IMPACT STATEMENT:** The accurate and precise SQA-iOw device offers a significant opportunity for point-of-care fertility practitioners to provide immediate, on-site, and reliable comprehensive semen analysis to their patients.

#### First Presenting Author

#### ***Presenting Author***

---

Natan Bar-Chama, M.D.

**Email:** nbarchama@gmail.com -- Will not be published

Reproductive Medicine Associates of New York Department of Urology The Mount Sinai  
Hospital  
635 Madison Ave  
Fl 10  
New York NY 10022-1009  
USA

Within the past 2 years, have you or your spouse/partner had any potential COI?  
Yes

Organization Name	Relationship Type	Who has this Relationship?
genentech/ roche Pharmaceuticals	Grant Recipient Relationship Began - Sunday, August 13, 2023 Relationship Ended - Honoraria Relationship Began - Sunday, June 11, 2023 Relationship Ended -	Self
WINFertility	Paid Consultant Relationship Began - Sunday, January 1, 2023 Relationship Ended -	Self

Signature: Natan Bar-Chama



Second Author

Lev Rabinovitch, PhD  
**Email:** levr@mes-ltd.com -- Will not be published

Medical Electronic Systems  
 6345 Balboa Blvd #185  
 Encino CA 91316  
 USA

Within the past 2 years, have you or your spouse/partner had any potential COI?  
No  
 Signature: Lev Rabinovitch, PhD

**CV Upload:**

Lev Rabinovitch Resume.pdf

---

**If necessary, you can make changes to your abstract submission until Thursday, April 25, 2024 at 5:00 pm (EDT).**

To access your submission in the future, use the link to your user portal from one of the automatic confirmation emails that were sent to you during the submission.

Or point your browser to <https://asrm.confex.com/asrm/2024/gateway.cgi>

You will be prompted to login with your ASRM account prior to accessing the user portal. If you do not yet have an ASRM account, the screen will redirect you to the site where you can register for a new account.

When registering for a new ASRM account:

- Please use the email address that is associated with the submission and your first and last name as they appear on the submission when creating this account.
- An ASRM account for login must be completed before you can access the user portal.

Any changes that you make will be reflected instantly in what is seen by the reviewers. You DO NOT need to go through all of the submission steps in order to change one thing. If you want to change the title, for example, just click "Title" in the abstract control panel and submit the new title.

When you have completed your submission, you may close this browser window.

If you would like to submit another abstract, click [here](#).

[Tell us what you think of the abstract submission process](#)

[Home Page](#)