

## SQA-iO Test Parameters

WHO 5 <sup>th</sup>	WHO 6 <sup>th</sup>
CONCENTRATION (M/ml)	CONCENTRATION (M/ml)
TOTAL MOTILE PR + NP (%)	MOTILITY (%)
PROGRESSIVE PR (%)	PROGRESSIVE (%)
NON-PROGRESSIVE NP (%)	RAPIDLY PROGRESSIVE (%)
IMMOTILE IM (%)	SLOWLY PROGRESSIVE (%)
MOTILE SPERM CONC.* (M/ml)	NON-PROGRESSIVE (%)
PROG. MOTILE SPERM CONC.* (M/ml)	IMMOTILE (%)
NORMAL FORMS (%)	NORMAL FORMS (%)
SPERM MOTILITY INDEX*	MOTILE SPERM CONC.* (M/ml)
SPERM # (M/ejac)	PROG. MOTILE SPERM CONC.* (M/ml)
MOTILE SPERM* (M/ejac)	RAPID PR. MOTILE SPERM CONC.* (M/ml)
	SLOW PR. MOTILE SPERM CONC.* (M/ml)
	FUNCTIONAL SPERM CONC.* (M/ml)
	VELOCITY (VCL)* (mic/sec)
	SPERM MOTILITY INDEX*
	SPERM # (M/ejac)
	MOTILE SPERM* (M/ejac)
	PROG. MOTILE SPERM* (M/ejac)
	FUNCTIONAL SPERM* (M/ejac)
	MORPH NORMAL SPERM* (M/ejac)