



# Certificate of Analysis

Product Name	<b>QwikCheck Liquefaction</b>
Product Number	0900
Product Lot	230709_0324
Storage Temp	Room temperature

<b>Test</b>	<b>SPECIFICATION</b>	<b>RESULTS</b>
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	> = 1200	1366
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: March 2025

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature:         S.Z          
                    Production Lab Lead Tech

Date:         27<sup>th</sup> March 2024