

QwikCheck™ Beads Batch Release Form

BATCH #	031223	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	December-24	LOT#	031223001	LOT#	031223002	LOT#	031223003		
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated		48	+/-	6.7	25	+/-	5.0	0.0	0.0
SQA-V Visualization		46	+/-	11.5	26	+/-	6.5	0.0	NA
Hemocytometer		49	+/-	12.3	25	+/-	6.3	0.0	NA
MAKLER		49	+/-	12.3	24	+/-	6.0	0.0	NA
Fixed Coverslip		50	+/-	12.5	26	+/-	6.5	0.0	NA

S.Z, Production Lab Lead Teach

S.Z

18-Dec-23

Signature

Date

M.D, CEO

Marcia Deutsch

18-Dec-23

Signature

Date