



## Certificate of Analysis

|                |                               |
|----------------|-------------------------------|
| Product Name   | <b>QwikCheck Liquefaction</b> |
| Product Number | 0900                          |
| Product Lot    | 219634_0326                   |
| Storage Temp   | Room temperature              |

| Test               | SPECIFICATION      | RESULTS |
|--------------------|--------------------|---------|
| Appearance (Color) | White to off-white | White   |
| Appearance (Form)  | Powder             | Powder  |
| Units/mg protein   | $\geq 1200$        | 1252    |

|                 |                                                                                                              |
|-----------------|--------------------------------------------------------------------------------------------------------------|
| Unit Definition | That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition |
|-----------------|--------------------------------------------------------------------------------------------------------------|

Expiration Date: MAR/2027

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature:   
Clinical Director

Date: 02-MAR-2026