



Technical Release Bulletin

QwikCheck GOLD Sperm Quality Analyzer (Software version 1.00)

Issue date: May 15, 2011

Subject: Differences between the QwikCheck GOLD and SQA-V

Intended use:

The QwikCheck™ GOLD Sperm Quality Analyzer is a standalone system for analyzing fresh, washed and frozen human semen samples.

Hardware Differences vs. SQA-V:

- No visualization compartment for viewing semen samples.
- No connection to the PC and no V-Sperm software.
- No built-in printer – an external label printer is available.

Software Differences vs. SQA-V:

- There is no option of testing post-vasectomy samples as they require a visualization system.
- No archive for semen and control test results.
- Only one concentration standard: Hemocytometer.

Compliance Date: Effective May 15, 2011

Authority: Lev Rabinovich, PhD CTO

Distribution: All distributors except U.S. on demand



www.mes-global.com

QwikCheck™GOLD vs SQAV Differences_15_MAY_2011