

TECHNICAL BULLETIN: SQA-V GOLD ACCESSORY STORAGE CONDITIONS

Thursday, August 21, 2014

BACKGROUND:

The storage conditions for the SQA-V GOLD accessories and kits are provided with each product separately when/if controlled conditions are required and are specified in the Material Safety Data Sheet (MSDS) for each product. This technical bulletin provides an overall list of the SQA-V GOLD accessories storage conditions which are provided by the manufacturer.

OVERVIEW:

Product Name	Storage Conditions
SQA-V GOLD System	<ul style="list-style-type: none">Storage Conditions are: 15-38 °C, up to 80% humidity @ 31 °C.Operational Conditions are: Room Temperature*
SQA-V GOLD Testing Capillaries	<ul style="list-style-type: none">Storage Conditions are: Room Temperature*
SQA-V GOLD Cleaning Kit	<ul style="list-style-type: none">Storage Conditions are: Room Temperature*
QwikCheck™ Test Strips	<ul style="list-style-type: none">Storage Conditions are: 15-30°C (60-86°F) and humidity <30%, away from direct sunlight. Keep the bottle tightly capped. Avoid touching test pads.
QwikCheck™ beads	<ul style="list-style-type: none">Storage Conditions are: Room Temperature*
QwikCheck™ Liquefaction kit	<ul style="list-style-type: none">Storage Conditions are: Room Temperature*
QwikCheck™ Dilution kit	<ul style="list-style-type: none">Storage Conditions are: Room Temperature*After opening bottle: Refrigerate at +4 °C (~39.2 °F)
SQA-V GOLD Validation Kit	<ul style="list-style-type: none">Storage Conditions are: Refrigerate at +4 °C (~39.2 °F)

*Room temperature is a general term describing an indoor temperature range. MES follows the WHO guideline for this range as 20°C (68°F or 293K) to 25°C (77°F or 298K).

MANUFACTURER'S RECOMMENDATION:

It is recommended to strictly follow the storage conditions as specified in this technical bulletin and on the MSDS for each product provided by the manufacturer for all SQA-V accessories.

Dr. Lev Rabinovitch, CTO MEDICAL ELECTRONIC SYSTEMS

Distribution: All SQA Users



www.mes-global.com



www.a-tech-global.com