

QwikCHECK™ QC PIG

Sperm Quality Analyzer

QC a PORCINE AI sample
prior to insemination in less than 45 seconds!
Fast, simple and accurate!



- Target market: AI Insemination sites that want to QC samples
- Tests EXTENDED pig semen
- Prints results to label maker
- Extender settings
- Sample requirement: 500 µl
- Re-usable, washable testing capillary

EXTENDED SEMEN test results reported:

- Total Sperm Concentration M/ml
- % Motility
- # Sperm
- # Motile Sperm



“Remember, it all started with a Sperm!”