


**MEDICAL ELECTRONIC SYSTEMS
WWW.MES-GLOBAL.COM**

QwikCheck™ Beads Batch Release Form									
BATCH #	090921	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	September-22	LOT#	090921001	LOT#	090921002	LOT#	090921003		
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck gold Automated		47	+/-	6.6	26	+/-	5.2	0.0	0.0
SQA-V Visualization		47	+/-	11.8	26	+/-	6.5	0.0	NA
Hemocytometer		43	+/-	10.8	24	+/-	6.0	0.0	NA
MAKLER		42	+/-	10.5	26	+/-	6.5	0.0	NA
Fixed Coverslip		52	+/-	13.0	26	+/-	6.5	0.0	NA

Sapir Zrihan, Clinical research assistant



Signature

23_SEPT_21

Date

Frida Cohen, Clinical Manager



Signature

23_SEPT_21

Date

Marcia Deutsch, CEO



Signature

23_SEPT_21

Date