



QwikCheck™ Liquefaction Kit Batch Release Form

QwikCheck™ Liquefaction Kit produced: FEB 2022
Month Year

Chymotrypsin activity = 1200 units/mg enzyme

Notes:

- PASS criteria: Mean ≥ 1200 units/mg enzyme
- Certificate of Analysis (Faizyme) for each Lot Number should be attached to this form.

Product insert	100		
Vial labels	2000		
Box labels	100		
QC labels	100		
QC Approval date	01/Feb/22		
Preparation date	02/Feb/22		
Release date	15 th Feb 22		
Signature	S.Z		

* 0900; 20 vials.

Batch number: 210911_0222

Lot Number (Faizyme): 210911

Exp. Date: FEB 2023
Month Year

Clinical Research Assistant

15th Feb 2022

Date