

**QwikCheck™ Low Range Beads Batch Release Form**

BATCH #	090322	Level 1			Level 2			NEGATIVE CONTROL	
EXP DATE:	09-Dec-22	LOT#	090322001	LOT#	090322002	LOT#	090322003		
For use with:		Conc.			Conc.			Conc	MSC
SQA / QwikCheck™ gold Automated		15	+/-	6.0	3	+/-	2.0	0.0	0.0
SQA-V Visualization		15	+/-	6.0	3	+/-	2.0	0.0	NA
Fixed Coverslip		15	+/-	6.0	3	+/-	2.0	0.0	NA
Hemocytometer		14	+/-	5.6	NA	+/-	NA	0.0	NA

Victoria Goldenberg

Clinical research assistant



Signature

13<sup>th</sup> March 2022

Date

Frida Cohen, Clinical Manager



Signature

13<sup>th</sup> March 2022

Date

Marcia Deutsch, CEO

*Marcia Deutsch*

Signature

13<sup>th</sup> March 2022

Date