

QwikCheck™ QC, test and semen treatment kits

*QwikCheck Kits
are FDA and
CE certified*

QwikCheck™ Kits

FDA and CE certified QC, test and treatment kits can be used to liquefy and dilute semen samples, test for pH and leucocytes and provide a known target value for sperm concentration against which laboratories can test the SQA-V calibration, the proficiency of manual sperm counting procedures and the calibration of sperm counting chambers.

QwikCheck™ Beads QC Kit

Quality control latex beads for testing sperm concentration. For use on both manual and automated systems:

- SQA-V, SQA-Vb, SQA-Vp, SQA-Ve
- Neubauer hemacytometers
- Makler counting chambers
- Fixed cover slip chambers (like MicroCell™)

Each kit contains:

- 3 Levels of QC material: High, Low and Zero
- 5cc of each level of QC material
- Instructions for use on all types of chambers and automated systems

Item #0200

QwikCheck™ Liquefaction Kit

Used to quickly liquefy viscous semen samples.

Each QwikCheck™ Liquefaction kit contains:

- 20 single dose treatment vials
- Instructions for use

Item # 0900

QwikCheck™ Dilution Kit

Used to dilute semen samples without impacting motility.

Each QwikCheck™ Dilution kit contains:

- 50 CC sterile dilution media
- Instructions for use

Item # 0800

QwikCheck™ Test Strips

FDA approved semen WBC and pH reagent test strips.

Each QwikCheck™ Test Strip kit contains:

- 100 test strips
- Product insert

Item # 0700

