

QwikCHECK™ Gold Sperm Quality Analyzer

The QwikCHECK™ Gold sperm analyzer is an affordable medical device that performs a complete quantitative evaluation of semen quality in less than 2 minutes. This high performance system incorporates technology in electro-optics and computer algorithms to provide a quick, precise and accurate automated semen analysis.

*Automated
Semen Analysis in
less than 2 minutes!*

Features

- Runs fresh, frozen and washed human semen samples
- Disposable testing capillary
- Easy to use – menu driven testing
- CE certified
- Self-testing, self-calibrating
- Runs latex beads or stabilized sperm for QC

Automated Test Results

- Sperm Concentration
- % Motility
- % Normal Morphology
- Motile and Functional Sperm Concentration
- Progressively Motile Sperm Concentration
- Velocity
- Sperm Motility Index

Dynamic Range of the QwikCHECK™ Gold

SAMPLE	SPERM CONC.	MSC	Motility
FRESH	2-400 or < 2 M/ml	0.2-400 or <0.2 M/ml	0-100%
WASHED	2-200 or < 2 M/ml	0.2-200 or <0.2 M/ml	0-100%
FROZEN	Not reported	0.2-200 or <0.2 M/ml	Not reported

