## QwikCheck

## Advanced Testing Kits



## QwikCheck | DFIKit

In-vitro Diagnostic Kit for **DNA Fragmentation**Reporting via the SCD Method.



## DNA Fragmentation Index Kit

- Complete Kit Includes supplies for 10 Tests
- ✓ Saves slide prep time: <25 Minutes vs >50 Minutes
- ✓ WHO 6th Edition Compliant: Small, Medium, & Large Halos
- Affordable and easy to use with clear step-by-step instructions
- Compatible with the SQA-Vision Analyzer or any Microscope



Learn More

QwikCheck | DFI Kit is an in-vitro diagnostic tool for reporting sperm DNA Fragmentation Index (DFI) using the Sperm Chromatin Dispersion (SCD) method. It simplifies slide preparation for accurate sperm fragmentation assessment based on the 6th Edition (five categories) methodology. Halo sizes relative to sperm heads can be evaluated with a microscope or the SQA-Vision Automated Sperm Quality Analyzer. The kit includes reagents for 10 tests and detailed instructions.