

# SQA<sup>7</sup>

## Testing Capillaries Low Volume



Fertility Testing with  
10-Microliter Samples



*"Remember, it All Started  
... With A Sperm!"*

Scan QR to



Learn More

The **Low Volume Capillary (LVC)** is designed to test washed, frozen, and other types of fertility sample preparations using only 10 microliters of sample. Using our **SQA-Vision** or **SQA-iO** Analyzer, the testing capillary can assess motile sperm count, progressive motile sperm count, total motile sperm count, total progressive motile sperm count, and velocity, and much more.

For use with the **SQA-VISION<sup>7</sup>** & **SQA-iO<sup>7</sup>**