Certificate of Analysis

Product Name	**QwikCheck Liquefaction**
Product Number	0900
Product Lot	235992_1126
Storage Temp	Room temperature

Test	**SPECIFICATION**	**RESULTS**
Appearance (Color)	White to off-white	White
Appearance (Form)	Powder	Powder
Units/mg protein	>= 1200	1228
Unit Definition	That amount of enzyme causing a decrease in absorbance at 237nm of 0.0075 per minute at 25°C Unit Definition	

Expiration Date: NOV/2027

This document confirms that these parameters have been tested in a QC laboratory and have been satisfactorily completed.

Signature: _____ Date: 01-Nov-2026

Production Laboratory Assistant Manager

MES-0186-03